



## AP-450 Ceiling Lift Portable Lifting Solution for Any Room

## AP-450 Specifications

## Part Numbers

**Motors** 

Lift Motor: 24 VDC

Safe Working Load (SWL)

204 kg (450 lb)

Dimensions (L x W x H)

365 x 209 x 190 mm (14 x 8.2 x 7.5 in)

**Unit Weight** 

5.5 kg (12.2 lb)

**Lift Case** 

Flame Retardant ABS

Safety

ActiveSense Tilt Detection

SafeGuard Hook

Antimicrobial Strap with Wear Indicator Marks

- O

Emergency Stop

**Emergency Lowering** 

Low Battery Indicator

Standards

IEC 60601-1:2005/A1:2012

IEC 60601-1-6:2010/A1:2013

IEC 60601-1-11:2015

IEC 62366:2007/A1:2014

ISO 10535:2006

IEC 60601-1-2:2014 (Ed 4.0)

**Hand Control** 

Digital

IPx6 and IPx7 Rating

**Service Life** 

10 Years or 22,500 Cycles

**Rated Performance** 

80 lifts at Max Load

**Maximum Sound Level** 

50 dB

**Batteries** 

Single 24V 4.0AH

(20 cell 1.2V) NiMH battery

**Charger Model** 

PS2403SRM26

**Charger Input** 

100-240 VAC 0.5 Amps, 50-60 Hz

**Charge Output** 

24 VDC, 1.33 Amps

**Lifting Range** 

Up to 1960 mm (77 in)

**Duty Cycle** 

1 Min "On"

9 Mins "Off"



AP-450: Lift Up to 204 kg (450 lb) 303060: Portable Ceiling Lift

We're Trusted.

Over 275,000 ceiling lift installations across North America

We're Experts.

Handicare handles the installation and provides **FREE on-site training** for your staff.



Take your safe patient handling to the next level - or anywhere - with portable ceiling lift technology.

