




**EASYTRACK 4 POSTS WITH
MOVING 3RD TRACK**
**Pressure fit system for portable ceiling
lifts**

ARJOHUNTLEIGH
GETINGE GROUP



...with people in mind

The Easytrack is a pressure fit system requiring no permanent installation.

The 4 posts configuration with moving 3rd track of the Easytrack system works on two X-Y axes while optimizing lifting options anywhere within the transfer area.

	Part #	Description
1	9440030	Easytrack - 4 Post (Ceiling lift sold separately)
2	700.10250	X-Y Tracks 3000 mm (118")
3	700.11005	KWIKtrak autolock track stopper
4	700.10650	Trolley for portable lifts/Easytrack
5	200.11050	Endcap
6	700-11225	X-Y mobile trolley
7	200.10730	Drilling JIG for track hook
8	700.10600	Clip-on level
9	203.03280	Allen key 6mm
	203.02505	Allen key 4mm
10	200.10720	Round endcap
11	700.10400	Top plate for post
12	700.10350	Round foot plate for post
13	700.10230	Post strip (2 pieces)
14	700.10500	Standard post
	700-08310	Reacher - 610 mm (24")
15	700-08320	Reacher - 910 mm (36")
	002.10450	Easytrack transport bag (Custom bags are required if tracks are greater than 1830 mm (72 in))



Product Information

Max. safe working load (SWL)	200 kg (440 lb)
Weight	57 kg (125 lb)
Range	
Length of rail	3000 mm (118")
Length of post	2130 - 2745 mm (84" to 108")
Material	Aluminium
Ceiling plate	510 mm (20")
Floor plate	120 mm (4.8")
Compliant with:	ISO 10535:2006

Product dimensions (#9440030)

