

# THE STATE OF THE ART IN MOBILITY



A person is seated in a white shower chair with armrests inside a shower stall. The person's hands are resting on the armrests. The shower stall walls are made of white tiles. The image has a light teal overlay.

**YOUR BRIDGE TO A MORE INDEPENDENT LIFE**

# WHY NUPRODX



## HISTORY

Nuprodx Mobility founder Bruce Hammer understands the need for mobility solutions as well as anyone. Bruce, a former professional motorcycle racer, was ranked second in the United States when his life took an unexpected turn in 1982. Bruce crashed during a race and became a quadriplegic after breaking his neck at the sixth/seventh cervical vertebra. Bruce, who also had a background in architecture, eventually began designing and creating products that would help him travel more independently. Inspired by the freedom that the products provided, Bruce founded Nuprodx in 1998 to offer solutions to others. Since then, Nuprodx has grown to produce the most complete, most innovative, and highest quality in-home and travel solutions on the market.

*\*US Patent 7,690,055 B2*

## MISSION

Nuprodx Mobility improves safety in and around the bathroom for seniors and people with disabilities, fostering a more independent life. With the most innovative design and made-in-the-USA manufacturing, Nuprodx Mobility delivers the highest quality and the most complete solutions for in-home use and travel.

## KEY BENEFITS

- Patented bridge technology\* reduces the number of transfers for safety and convenience.
- Eliminates need for expensive remodeling to make virtually any bathroom accessible.
- Made-to-order in the USA to meet just about any special requirement.
- Modular design allows for easy disassembly and portability for travel.



# BUILT TO ORDER IN THE USA



*Built to order to meet customer needs*



*Precision-machined for stability and durability*

Nuprodx Mobility designs and manufactures our products in our facility in Livermore, California, from the precision-machining of components to final assembly. Because every system is built to order, Nuprodx Mobility has the ability to choose from various modular parts, allowing us to configure systems that meet just about any special requirement. That means Nuprodx systems can fit most any toilet, tub and shower and adapt to any bathroom environment. In addition, Nuprodx chairs can be adjusted and customized to provide users of all shapes and sizes with amazing comfort and convenience.



*Made with pride in the USA*

# TRANSFERS MADE SIMPLE

Nuprodx Mobility gets your day started right. With just one bedside lateral transfer, Nuprodx users are able to move around the home freely, avoiding dangerous bathroom transfers onto toilets and into tubs.



*Just one bedside lateral transfer*



*Easy transport*



*Slides over bridge and into tub to eliminate transfer*



*Rolls over toilet to eliminate transfer*

# SOLUTIONS FOR ANY BATHROOM

“A \$3000 solution to a \$30,000 problem” - Cara Colgin-Smith, Southwestern Medical Reps



*Tub Slider system  
(MC6000Tilt shown)*



*Shower Slider System  
(MC6000RS shown)*

Nuprodx Mobility has a solution for just about any bathroom environment. With Nuprodx you can use your existing tub or shower without the need for an expensive and disruptive re-model. Tub Slider Systems can be customized to navigate virtually any tub shape and size. The patented bridge allows users to slide into the tub without a manual transfer. Shower sliders fit through narrow shower entries and over step-in shower thresholds. Barrier-free shower solutions include large-wheel, shower-ready chairs for independent users who self-propel.



*Barrier-free shower solutions  
(MC4020RX shown)*

# MEETING CUSTOMER NEEDS

## Seniors



*Extends independence of seniors, including in their own home (MC6000 shown)*

## Pediatric



*Easily modified to grow with children (MC6000TiltPed shown)*

Nuprodx Mobility meets our diverse customer needs. In addition to serving adults with disabilities, Nuprodx helps seniors remain in their own homes with more independence. The pediatric line is designed to grow with a child, with easy adjustments and modifications that can save thousands of dollars compared to investing in a new system. In all cases, Nuprodx designs its products with portability and travel in mind, giving additional meaning to the concept of mobility.

## Travel



*Easy to pack for travel and for peace of mind*

# TUB SLIDER SYSTEMS

Don't remodel it - bridge it



MC6000

Slider systems provide a commode chair, tub slider system, and in-home mobility chair (6000 series). Requiring just one bedside transfer, tub sliders eliminate the need for transfers onto toilets and into tubs in the dangerous bathroom environment. One transfer and you are ready to roll!

Tub Slider chairs slide easily over and across tub walls without a manual transfer. After the chair slides into the tub, the patented connector bridge can be removed, and the shower curtain can be closed for privacy.

## Key Benefits:

- Requires only one bedside transfer, providing safety and convenience.
- Modular design and built to order to fit most bathroom spaces.
- High quality construction for home reliability and key models that can be packed for travel.



# TUB SLIDER PRODUCT DETAILS



## MC5300: Bench Slider

Integrated design allows for a more independent user to slide across bridge and into tub with little or no assistance. Intuitive design that can be set up once and left in place. Perfect solution for seniors.



## MC6000: Best Seller

Our best seller. Commode chair, tub access slider system and in-home mobility chair. Locking cantilever arms swing up to ease transfers.



## MC6000TX: Road Warrior

All the benefits of the MC6000 and packs into the supplied roller suitcase. Proven road warrior for travel.



## MC6200: Extra Support

All the benefits of the MC6000 with a larger 20"x20" seat frame and a 400lb weight capacity.



## MC6000Tilt: Comfort Tilt

A center-of-mass, tilt-in-space, full-featured model with head support for a person with a higher-level disability, providing maximum comfort and stability.



# SHOWER SLIDER SYSTEMS

Slide in safe



MC6000RS

The perfect solution for step-in showers, including those with narrow door entries. Commode chair, in-home mobility chair, and shower access slider system that eliminates the need for dangerous transfers in the bathroom. Shower Slider chair enters backward through shower door to maintain a slim profile, slides over and across step-in-shower lip using patented connector bridge, then swivels to face shower head and locks in place.

## Key Benefits:

- Requires only one bedside transfer, providing safety and convenience.
- Slim profile and reverse orientation allow entry through narrow shower doors.
- Shower turntable allows an easy rotation to face shower head and controls.



# SHOWER SLIDER PRODUCT DETAILS



## **MC6000RS: Shower Ready**

With locking cantilever arms that swing up to ease transfers and innovative design to enter narrow step-in showers.



## **MC6000RSTilt: Comfort Tilt**

All the benefits of the 6000RS with tilt-in-space capability for maximum support and comfort.

# BARRIER-FREE SHOWER SOLUTIONS

Roll in, roll out



MC4024

Designed for those with barrier-free showers, the 4000 series models require just one bedside lateral transfer, eliminating the need for transfers onto toilets and into tubs in the dangerous bathroom environment.

Caster-wheel-only models are designed to enter the shower with caregiver assistance, while models with larger wheels are perfect for more independent users who can self-propel.

## Key Benefits:

- Requires only one bedside transfer, providing safety and convenience.
- Variety of wheeled solutions for smooth rolling on tile, carpet, and other home surfaces.
- Self-propelled line includes narrow footprint models for tighter home environments.

## MC4020RX: Front Mounted:

For more independent users. Mid-size 20" wheels with over-sized hand rims provide more narrow profile for navigating tight spaces in the home. Wheels are front-mounted to provide back clearance for larger toilet tanks.



# BARRIER-FREE SHOWER PRODUCT DETAILS



## **MC4000: Easy Roller**

For those with caregiver assistance, this chair can be rolled right into a barrier-free shower. With locking cantilever arms that swing up to ease transfer and with caster wheels that can be locked in place.



## **MC4000TX: World Traveler**

All the benefits of the MC4000 and the ability to be easily disassembled and packed into a supplied FAA compliant carry-on travel bag.



## **MC4000Tilt: Comfort Tilt**

All the benefits of the MC4000 with tilt-in-space capability for maximum support and comfort.

# STATIONARY CHAIRS

Easy and secure transfer



*MC3000TX with included transfer cushion*

## Stationary Shower & Commode Chairs

Designed for those with good mobility and developed transfer skills who want a more portable solution for commodes and showers or tubs.

### Key Benefits:

- Lineup offers lightweight and portable models such as the MC3000TX that can be packed into a convenient lap-sized carry-on travel bag.
- Innovative design with clamping system secures chair legs firmly in place, for rigidity and solid feel you would expect to find in a heavier, non-portable welded chair.
- Strong and durable chairs with high standard of weight-bearing, corrosion prevention, and resistance to water absorption.

### MC3000TX [left]

The included transfer bench that rests on top of the tub wall and the swing-away and removable arms allow user easy transfer from wheelchair onto bench and into shower chair.



### MC3000TX: Frequent Flier

All the quality of the MC3000 plus a handy transfer bench. Cloth back for even better portability— allows it to be packed down into supplied lap-size carry-on travel bag (18" x 9.5" x 9.5").



# STATIONARY CHAIRS PRODUCT DETAILS



## **MC3000: Base Product**

Our base product with all the quality you expect from Nuprod Mobility— a lightweight, durable and highly functional shower/commode chair that can go nearly anywhere.



## **MC3000Tilt: Comfort Tilt**

Quality of the MC3000 plus tilt-in-space capability for maximum support and comfort.



## **MC3200: Extra Support**

Same benefits of the MC3000 plus larger seat frame and 400lb weight capacity for customers needing more space and strong support.

# PEDIATRIC CHAIRS

Grow with your child



*MC6000TiltPed*

With its modular construction, the pediatric line is designed to grow with your child. With minor adjustments and/or updates to the seat and backrest cushions, you can enlarge the chair as the user grows, saving thousands of dollars compared to investing in a whole new system.

## **Key Benefits:**

- All the benefits of the adult systems, including durability, stability, and fewer and safer transfers.
- Adjustable seat width and depth to grow with a child.
- Modular design allowing for updates to seat and backrest cushions for more significant growth and changes in size.

## **MC6000TiltPed [left]**

Tub slider system requires just one bedside lateral transfer, eliminating the need for transfers onto toilets and into tubs in the dangerous bathroom environment. Variable seating components let it grow with the user.



# PEDIATRIC CHAIRS PRODUCT DETAILS



## **MC6000TiltPed: Fully Loaded**

A center-of-mass, tilt-in-space, full-featured model with head support, providing maximum comfort and stability.



## **MC6000Leckey: Mesh Comfort**

Tub-slider system with unique mesh chair that provides the ultimate comfort for showering.



## **MC4000TiltPed: Roll-in Comfort**

Designed to be rolled into barrier-free showers with caregiver assistance, a center-of-mass, tilt-in-space, model with head support for maximum comfort.

# ACCESSORIES



*Drop rail aids access for improved hygiene*



*Tub-side transfer bench eases bath entry*



*Splash guard*



*Supportive calf panel*



*Custom cushions*



*Shower legs*



*Commode pan (standard mount)*



*Optional bridges*



*Solid seat cushions available*



*Floor-support transfer bench*



*Commode pan with lowering brackets*



*Side hip pads*



*Optional seat belts*



*Suction cup feet and crutch tip feet options*



*Arm Troughs*



*Full selection of carrying cases*



THE STATE OF THE ART IN MOBILITY



©2019, NuprodX Inc.

Any or all specifications subject  
to change without notice.

