



MULTICHAIR 3000Tx
(carrying case included)


The lightest, smallest, and best portable tub transfer bench / commode chair available. So small that when folded and packed it can fit on your lap, and with its standard-equipment carrying case, it complies with FAA requirements for carry-on luggage. It's so compact, it can fit underneath most standard wheelchairs.

Features

- Includes swing-away transfer cushion kit
- Cantilever fold-back / locking padded chair arms
- Self-skinning padded seat cushion with 1/2" waterfall-style foam ring around the commode opening and slot help protect the ischial and coccyx areas to reduce risk of pressure sores
- Waterproof, fast-drying nylon cloth seatback
- Black nylon carrying case included (9" x 9.5" x 18")
- Capable of going through most airport X-ray machines
- Won't rust or corrode—constructed of aircraft-grade aluminum, brass, stainless steel and durable plastics

Specifications

- Approximate weight: 14 lbs.
- Weight capacity: 300 lbs.
- Passed ANSI / BIFMA X5.1 test performed by Gaynes Labs, Inc. (a national testing laboratory). Weights of 250 lbs. and 300 lbs. were dropped on the chair, with no structural breakage or loss of serviceability



3000Tx

18.5" overall depth

* Illustration shows available A3 armrests.