



Hillrom™

Advanta™ 2 Bed
Rental offering

SIMPLY INTUITIVE

We listened to nurses from Med-Surg units everywhere to integrate a number of features that facilitate an intuitive, empowering experience.

At Hillrom, we understand that Med-Surg caregivers must possess a multitude of capabilities to treat patients effectively. We made sure the Advanta™ 2 bed can help you navigate your ever increasing daily challenges.

The results: the Advanta 2 bed integrates several features that facilitate an intuitive, empowering experience for caregivers and patients alike.

INTUITIVE CAREGIVER EXPERIENCE

- Four-corner brake and steer for easy mobility
- One-button Dining Chair® feature places patients at a low height
- Built-in bed extender lets a single caregiver add length to the bed
- Integrated Line Manager keeps IV lines untangled and out of the way

PROTECT PATIENTS FROM SKIN BREAKDOWN

Designed by experienced clinical professionals, our comprehensive Safe Skin® program for the prevention and treatment of skin breakdown helps improve outcomes with the goal of making them repeatable and sustainable.

The bed's open architecture also lets you stretch your budget by choosing virtually any available support or therapy surface. By empowering patients with features that enhance their autonomy and facilitate their independence, the Advanta 2 bed improves their experience and is designed to deliver comfort.

EMPOWERING PATIENT EXPERIENCE

- Care Grip patient egress assist handle
- In-rail storage compartment
- Siderail patient positioning controls are freely accessible to the patient

A healthcare professional, a woman with brown hair in a ponytail wearing light blue scrubs and a stethoscope, is assisting a male patient. The patient is wearing blue hospital scrubs, glasses, and white socks. He is sitting on the edge of a white hospital bed, leaning forward with his right arm resting on the bed's control panel. The nurse is standing to his left, holding his right arm. The bed has a control panel with various buttons and a digital display. The background shows a hospital room with wooden cabinets and a tiled floor.


The Care Grip patient egress assist handle is designed to facilitate safer bed entry and exit by providing ergonomically correct leverage for patients.

SIMPLY RELIABLE

Patient safety is the number one concern of every caregiver. That's why we created these state-of-the-art features to help the caregiver avoid potential safety issues and support fall prevention programs or protocols. The Advanta™ 2 bed is designed to promote safety at the bedside and throughout the facility.



One-button Dining Chair® position enables caregivers to place patient in the chair position while the bed remains at a low height.



SlideGuard® bed articulation technology helps support patient positioning.



The ergonomically correct Care Grip handle provides just the right grip for safer bed exit.



Point-of-Care® siderail controls are conveniently located in a central location, on both sides of the bed, where patient and caregiver interaction occurs.



The Electric Vascular Foot allows the caregiver to raise patient's feet without having to lift the weight of the patient's legs.



Convenient siderail storage provides the patient access to their cell phone, glasses or other personal items.

SIMPLY OUTSTANDING

The Advanta™ 2 Med-Surg bed is designed to help you best meet your daily challenges. Its intuitive and integrated functionality — like positioning, transporting, communication, monitoring, and alerts — deliver the features to help address your top priorities and enhance your experience as a caregiver.



Three-Mode Bed Exit System with three different sensitivities — Patient Movement, Bed Exiting, and Out-of-Bed — may help reduce false alarms with normal patient movement.



SafeView® Alert System keeps staff informed of a patient's bed status.



The Advanta 2 bed connects to send data to NaviCare® Nurse Call, EMRs or other clinical applications.



SIMPLY EMPOWERING

9



7

8



- 1 Line-of-Site® Feature
 - 2 Chair Position
 - 3 Care Grip Handle
 - 4 Three-Mode Bed Exit System (Patient Position Monitor)
 - 5 Integrated Scale
 - 6 OneStep® Operation Siderails
 - 7 Point-of-Care® Siderail Controls
 - 8 Electric Vascular Foot
 - 9 Emergency CPR
- Brake Alarms
- Battery Backup
- Built-in Bed Extender
- SlideGuard® Technology Frame
- Ability to send data to Nurse Call, EMR or other clinical applications

TECHNICAL SPECIFICATIONS

Frame Width	40" (101.6 cm)
Frame Length with Head/Footboard (w/Bumpers) Normal Extended	91" (231 cm) 94" (238.9 cm)
Sleep Deck Width Normal Length Extended Length	36" (91.44 cm) 84" (213.4 cm) 87" (221 cm)
Bed Clearance Under Frame	5.00" (12.70 cm)
Scale System Accuracy Weigh In Any Position Scale Repeatability	+/- 2.2 lbs or 1.0% of patient weight w/e is greater yes 0.2% of patient weight
Low Position Floor to Deck	15.75" (42.55 cm)
High Position Floor to Deck	32.5" (82.55 cm)
Max Head Elevation	65°
Max Trend and Reverse Trend	16°
Caster Size	6" (15 cm)
Safe Working Load	500 lbs



Hillrom™

**Healthcare can be unpredictable.
Your rental partner shouldn't be.**

For more information, please contact your Hillrom sales representative at 1-800-445-3730 or order online at rental.hillrom.com.

hillrom.com

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

© 2021 Hill-Rom Services, Inc. ALL RIGHTS RESERVED. 161274 rev 4 17-MAR-2021 ENG – US