

Geo-Mattress® 350


High-Value Option Meets Multiple Needs

Segmented top surface, 6" height and fluid-proof cover make the new 350 a great choice for heavier home care patients. Or, an economical yet effective option for long term care residents.

Item #	Description
67535-29	Geo-Mattress 350, 75"L x 35"W x 6"H
68035-29	Geo-Mattress 350, 80"L x 35"W x 6"H
68435-29	Geo-Mattress 350, 84"L x 35"W x 6"H
Medicare Approved: E0184 (Group I)	
Weight limit:	350 lbs. Mattress Weight: 21 lbs.
Warranty:	5 years, (2 years prorated) [applies in all care settings]

Description

6" Therapeutic Foam Mattress



All models are fire barrier version (FBV) as indicated by "-29" suffix, denoting certified conformance to NFPA 101 (Life Safety Code); Cal Tech Bulletin #117; 16 CFR Parts 1632 and 1633; and BFD IX-11.