

# Geo-Mattress® HC

## Your Best Bet in Group I Surfaces

Specifically designed to aid in pressure ulcer prevention in the home care setting. Geo-Mattress HC is cut from a single layer of high-quality, combustion modified foam in a fluid-proof, nylon and vinyl cover. Cross-cut top surface is segmented in three distinct zones. Folded and boxed to save freight and storage.

Note: This product does not include Heel Slope® feature.

Item #	Description
SP812-29	Geo-Mattress® HC, 75"L x 35"W x 5"H
SP810-29	Geo-Mattress® HC, 80"L x 35"W x 5"H
SP811-29	Geo-Mattress® HC, 84"L x 35"W x 5"H

**Medicare Approved:** E0184 (Group I)

**Weight limit:** 250 lbs.      **Mattress Weight:** 17 lbs.

**Warranty:** 2 years, non-prorated [applies in Home Care use only]

## Description

Therapeutic Home Care Foam Mattress



**All models are fire barrier version (FBV)** as indicated by "-29" suffix, denoting certified conformance to NFPA 101 (Life Safety Code); Cal Tech Bulletin #117; 16 CFR Parts 1632 and 1633; and BFD IX-11.