

Geo-Mattress® HC

Category: Therapeutic Home Care Foam Mattress

Specifically designed to aid in pressure ulcer prevention in the home care setting. Geo-Mattress HC is cut from a single layer of high-quality, combustion modified foam in a fluid-proof, nylon and vinyl cover. Cross-cut top surface is segmented in three distinct zones. Folded and boxed to save freight and storage.


Note: This product does not include Heel Slope® feature.



| Item # | Description |
|----------|-------------------------------------|
| SP812-29 | Geo-Mattress® HC, 75"L x 35"W x 5"H |
| SP810-29 | Geo-Mattress® HC, 80"L x 35"W x 5"H |
| SP811-29 | Geo-Mattress® HC, 84"L x 35"W x 5"H |

| | |
|--------------------|--|
| Medicare Approved: | E0184 (Group I) |
| Weight limit: | 250 lbs. |
| Warranty: | 2 years, non-prorated [applies in Home Care use only] |
| Mattress Weight: | 17 lbs. |

All models are fire barrier version (FBV) as indicated by "-29" suffix, denoting certified conformance to NFPA 101 (Life Safety Code); Cal Tech Bulletin #117; 16 CFR Parts 1632 and 1633; and BFD IX-11.



Geo-Mattress® 350

Category: 6" Therapeutic Foam Mattress

High-value option meets multiple needs. Segmented top surface, 6" height and fluid-proof cover make the new 350 a great choice for heavier home care patients. Or, an economical yet effective option for long term care residents.

| Item # | Description |
|----------|-------------------------------------|
| 67535-29 | Geo-Mattress 350, 75"L x 35"W x 6"H |
| 68035-29 | Geo-Mattress 350, 80"L x 35"W x 6"H |
| 68435-29 | Geo-Mattress 350, 84"L x 35"W x 6"H |

| | |
|--------------------|---|
| Medicare Approved: | E0184 (Group I) |
| Weight limit: | 350 lbs. |
| Warranty: | 5 years, (2 years prorated) [applies in all care settings] |
| Mattress Weight: | 21 lbs. |

All models are fire barrier version (FBV) as indicated by "-29" suffix, denoting certified conformance to NFPA 101 (Life Safety Code); Cal Tech Bulletin #117; 16 CFR Parts 1632 and 1633; and BFD IX-11.

Gel 80™

Category: Gel/Foam Mattress Overlay

The Gel-80™, an advanced Group 1 support surface gel/foam overlay, delivers effective pressure redistribution by encapsulating multiple gel bladders within a high-density sculpted support base. Features include a wipe-clean, fluid-impervious fabric top surface, elastic corner straps and one-piece construction. Folded and boxed for easy shipping and storage.



| Item # | Description |
|--------|---|
| GEL-80 | Gel/foam overlay for hospital-sized twin (80x35) mattress |

| | |
|--------------------|-----------------|
| Medicare Approved: | E0185 |
| Weight limit: | 300 lbs. |
| Dimensions: | 78"Lx34"Wx3.5"H |
| Warranty: | 1 year |
| Overlay Weight: | 46 lbs. |

Components conform to applicable California Technical Bulletin 117 standards.

