



*The perfect balance between
comfort and function*

Harmony

All Purpose Healthcare Reclining Chair

Standard Features

- A highly comfortable recliner offering excellent support and positioning built on an attractive yet durable tubular steel frame.
- Two 4" total lock swivel casters (diagonally positioned for secure locking), one front directional lock swivel caster for steering control during transport, and one rear swivel caster.
- Comfortable removable/replaceable cushions using high quality high resiliency foams and upholstered in contract/healthcare grade leatherette (Silververtex vinyl with bacteria-fighting defense or Crypton fabric available).
- Wide padded armrests for comfort during prolonged sitting.
- Occupant operated recline with adjustable back angle.
- High back design with push bar.
- Highly adjustable; features adjustable seat height, seat depth, initial seat angle, and backrest angle.
- Wide range of options and accessories permits configuring the chair for a variety of applications.



A lever locks the leg rest and recline angle in the desired position.

New Features

- Lever activated Locking Mechanism (C-C15 option) maintains the desired leg rest/recline position to optimize comfort and proper positioning.
- New graphite color silver nitrate impregnated powdercoat paint for improved infection control and outstanding durability.
- Silver nitrate impregnated Silververtex vinyl for improved infection control (optional).
- Now available in two widths, 17" (43cm) and standard 20" (51cm).



Chair shown with optional C-C15 leg rest/recline position lock.

Harmony Chair – Options and Accessories

More images showing different configurations can be viewed on our Website: <http://www.lpamedical.com>



C-C1



C-C2



C-C3



C-C4



C-C6



C-C7



C-C8



C-C9



C-C10



C-C11



AB-08H



AE-02



RA-03



RA-03A



RA-19



RA-24



RB-18

Item	Description
C-C1	Left and/or Right Lateral Tray
C-C2	Adjustable Elbow Supports
C-C3	Adjustable I/V Pole, Rear Installation
C-C3A	Adjustable I/V Pole, Front Installation
C-C3B	Adaptor for IV Pole with 1/2" Round Tube
C-C4	Elevating Foot Cushion
C-C5	Four (4) Swivel Wheels with Brake
C-C6	Self-Storing and Removable Footrests
C-C7	Side Comfort Cushions
C-C8	Opaque Activity Tray
C-C8-2	Activity Tray with Clear Window
C-C9	Lateral Support Bracket for Storing Activity Tray
C-C10	Pre-Cure Cushion Allows the Insertion of Specialty Pressure Management Gels/Foams
C-C11	Inclined Armrest for Blood Sampling
C-C15	Leg Rest/Recline Position Lock Mechanism
AB-08H	Contoured Headrest
AC-04	Leatherette Padded Tray Cover
AE-02	Catheter Drainage Bag Hook
RA-03	Lateral Safety Straps for Activity Tray
RA-03A	Full Table Security Strap - Secures Table Around Back
RA-19	Universal Perineal Abductor Kit with Black Fixing Mat
RA-23	4" (10cm) Lateral Trunk Supports with Backrest Canvas
RA-24	6" (15cm) Lateral Trunk Supports with Backrest Canvas
RB-18	Cervical Support

Technical Specifications

Model	C-207-17	C-207-20
Overall Width	25 1/2" / 65cm	28 1/2" / 72cm
Closed Length	36 1/2" / 93cm	
Opened Length	70" / 178cm	
Closed Height	46 1/2" / 118cm	
Seat / Floor Height	Adjustable from 17 1/2" to 20 1/2" (44 to 52cm)	
Seat Width	17 1/4" / 44cm	20 1/4" / 51cm
Seat Depth	Adjustable from 17" to 18 1/2" (43 to 47cm)	
Seat / Backrest Angle	Adjustable from 95° to 108°	
Weight w/o Accessories	75 lbs / 34kg	80 lbs / 36kg
Weight Capacity	250 lbs / 114kg	

Please note that technical specifications are subject to change without notice.

Upholstery – Contract grade leatherette or optional Crypton. Meets CA Technical Bulletin 117-2013.

The colors shown can be different from the real product.



Applications

The Harmony is the most versatile chair in its class thanks to its unique design. It offers excellent value for the money. With its exceptional comfort it can be used for blood sampling, specific treatments, general or geriatric care and prolonged sitting. Ideal for palliative care where comfort and position change is important. Its new optional leg rest/recline position lock sets it apart from common care recliners on the market. Furthermore, its design and new graphite powder coat make it more attractive and less clinical looking.

Our products – The complete line



Comp-Position
Positioning Chair



Dyn-Ergo
Scoot Chair



Evolution
Mobility Chair



Harmony
Care Recliner



Nomad
Mobile Comfort



Optimum
Posturo-Pedic
Positioning Chair



Wood
Thera-Glide®
Safety Glider



Decorative Metal
Thera-Glide®
Safety Glider



Transitional Style
Thera-Glide®
Safety Glider



Utility
Stacking Chair



Versatile
Adjustable Patient
Room Chair



Pre-Cure
Pressure Relief-Insert
Ready Cushion



Anti-Slip Surface
Posturo-Pedic
Positioning Aid

**Contact us for
more details
on our range of
products, or visit
our website.**