



**Sirocco  
Shower stretcher**



6100 2400

# Sirocco Shower stretcher

6100 2400

## Sirocco Shower stretcher

wall-mounted outlet left

- Maximum lifting capacity 150 kg.
- Electrically operated height adjustability.
- Easy-to-use hand switch.
- Suitable for bathrooms with minimal ceiling height of 230 cm.
- Incorporates wall-mounted battery charger.
- Battery easy to change.
- Emergency stop.
- Emergency lowering device.
- Length: 199 cm, width: 75 cm.
- Lowest lying position: 40 cm.
- Stretcher, head and foot boards made from polycarbonate.
- Mattress has 20 cm sill.
- Height adjustability mechanism can be operated by handcontrol.
- 2 hinged safety rails, coated stainless steel.

- The stretcher is constructed from coated steel with protective plastic trim.
- Ideal as shower stretcher and changing table.
- The height adjustability mechanism is powered by a 24 V actuator, fitted with a removable battery pack which can be wall-mounted outside the bathroom.
- The stretcher is fitted with a 1.75 m water drainage outlet hose fitted on right-hand side (standard) when viewed from the operating side. Longer outlet hoses available on request.
- Standard NEN 1010 to use safe in wet area.
- **Optional:** extra battery pack.

6100 2401

## Sirocco Shower stretcher

wall-mounted outlet right

