

ARJOHUNTLEIGH

GETINGE GROUP

CAREVO SHOWER TROLLEY



...with people in mind

Change the way you shower your patient

Carevo™ Shower Trolley from ArjoHuntleigh is an innovative hygiene solution for immobile patients.


The new design and intuitive interface provide an easy-to-maneuver and easy-to-use hygiene solution for both caregivers and immobile patients.






1 Integrated Control Panel

- Accessible at all times
- Simple, intuitive and efficient




2 Comfort Handles

- Enhanced sense of security
- Promotes participation and mobility
- Provides support



3 Flexizone

- Improved comfort
- Provides better pressure redistribution



4 Ergo Access Area

- Work closer to your patient
- Reduce static load
- Facilitates turning and repositioning



5 Elevated Back Support

- Optional
- Easier and efficient patient transfers
- Optimal patient comfort



6 Width Adjustment Option

- Greater comfort for a wide range of patient sizes
- Versatility to meet more of your showering needs



7 Integrated Shower Handle Holder

- Practical and intuitive design
- Convenient and easy-to-use for the caregiver



8 Optimized Head Pillow and softer mattress

- Superior patient comfort
- Durable, easy-to-clean and disinfect



9 Simple Daily Maintenance

- Easy mattress attachment and removal
- Easy-to-disinfect

Product Information	
Maximum Safe Working Load (SWL)	440 lbs (200 kg)
Maximum Patient Weight	401 lbs (182 kg)
Weight (<i>Carevo</i> with mattress)	181 lbs (82 kg)
Width (stretcher)	23" (586 mm)
Total Width (including side supports)	30" - 35" (759 mm - 888 mm)
Total Length	79" (2005 mm)
Length of Chassis (including wheel frame)	47 ½" (1203 mm)
Width of Chassis	27 ½" (700 mm)
Height of Chassis (over wheels)	8" (204 mm)
Height, high position from floor to stretcher	37 ½" (955 mm)
Height, low position from floor to stretcher	23 ¾" (605 mm)
Stroke Range	13 ¾" (350 mm)
Drain Hose Holders included	2 pieces
Maximum Back Support Angle (Optional)	15°
Equipped with head pillow and foldable side supports with safety lock	

ARJOHUNTLEIGH
GETINGE GROUP