

The solution for safer independent mobility, better positioning, improved comfort, enhanced dignity, and fall reduction.

Dyn-Ergo & Evolution

Standard Features

- Adjustable seat height.
- Adjustable seat depth.
- Gas spring assisted tilt-in-space with 8 mounting points makes tailored position change easier.
- Adjustable seat back angle.
- Adjustable armrest height.
- New graphite color silver nitrate impregnated powdercoat paint for improved infection control and outstanding durability.





Evolution S-250 Mobility Chair



Dyn-Ergo S-150 Scoot Chair

New Options

- 1 Fully adjustable elevating leg rest (optional).
- 2 Full-width padded footboard (optional).
- 3 Activity tray with clear window (optional).
- Contoured headrest (optional).
- Quick release seat depth adjustment (optional).
- Perineal abductor with universal black fixing mat (optional).
- Adjustable tension backrest (optional).
- Silver nitrate impregnated Silvertex vinyl for improved infection control (optional).







Dyn-Ergo & Evolution

Dyn-Ergo & Evolutions Chairs – Options and Accessories

More images showing different configurations can be viewed on our Website: http://www.lpamedical.com



Item	Description	S-150	S-250
SA-10	Swing-Away Self-Storing Removable Footrests	√	_
SA-10L	Left Footrest	√	_
SA-10R	Right Footrest	√	_
SA-15	Swing-Away Self-Storing Removable Footrests	_	√
SA-15L	Left Swing-Away Footrest	_	✓
SA-15R	Right Swing-Away Footrest	_	√
SA-20	Structural Lateral Supports	√	✓
SA-30	13" (33cm) Rear Wheels with Handbrakes	√	_
SA-35	Set of 4 Multi-Directional Total Lock 5" (13cm) Casters	√	_
SA-40	Activity Tray	✓	✓
SA-40-2	Activity Tray with Clear Window	√	✓
SA-45	Storage Support for Activity Tray	√	✓
SA-50	Cervical Support	√	✓
SA-60	Pre-Cure Cushion Allows the Insertion of Specialty Pressure Management Gels/Foams	✓	✓
SA-70	Heel Support Loops	√	√
SA-75	Height Adjustable I/V Pole	\checkmark	✓
SA-80	Increased Seat Depth Option Adjustable from $151/2"$ to $203/4" (39\text{cm}$ to $53\text{cm})$	✓	✓
SA-85	Quick Release Seat Depth Adjustment	√	√
SA-90	Increased Seat Height Option Adjustable from $161/2$ " to $211/2$ " (42cm to 55cm). Conflicts with SA-30 Option.	✓	_
SA-120	Full-Width Padded Footboard (Requires SA-10 Option)	√	√
SA-130R	Adjustable Elevating Leg Rest (Right)	√	√
	Adjustable Elevating Leg Rest (Left)	√	√
SA-215 (width)	Adjustable Tension Backrest	✓	✓
RA-03	Lateral Security Straps for Activity Tray	√	√
RA-03A	Full Table Security Strap for Activity Tray	√	√
RA-19	Perineal Abductor with Universal Black Fixing Mat	√	√
RA-23	4" (10cm) Adjustable Lateral Trunk Supports	√	√
RA-24	6" (15cm) Adjustable Lateral Trunk Supports	√	V
AB-08-H	Contoured Headrest	✓	✓
AE-02	Catheter Drainage Bag Hook	✓	√

Contract grade leatherette upholstery standard, Silvertex vinyl or Crypton fabric available. Conforms to CA Technical Bulletin 117-2013.

Technical Specifications

·	S-150 (Dyn-Ergo)					S-250 (Évolution)				
Model	16	18	20	22	24	16	18	20	22	24
Maximum Length			58"/147 cm					57"/145 cm		
Minimum Length			31"/79cm					32"/81 cm		
Overall Width	23"/58 cm	25 ½"/65 cm	27"/69cm	29"/74 cm	30 ½"/77 cm	25 ½"/65 cm	28"/71 cm	29½"/75cm	31½"/80cm	33"/84cm
Overall Width with SA-30 Option	27"/69 cm	29½"/75cm	31"/79cm	33"/84 cm	34½"/88cm	_	_	_	_	_
Maximum Height			54"/137 cm					541/2"/138 cm		
Minimum Height			36 ½"/93 cm					38"/97 cm		
Seat Width	151/2"/39cm	18"/46cm	19½"/50cm	21½"/55cm	23"/58cm	151/2"/39cm	18"/46 cm	19½"/50cm	21 ½"/55 cm	23"/58cm
Seat Depth	Adjustable from 13 ¾" to 18" (35 to 46 cm)									
Seat Height	Adjustable from 13 ½" to 18 ½"/34 to 47 cm				Adjustable from 16" to 21"/41 to 53 cm					
Armrest Height from Seat	Adjustable from 8" to 91/2"/20 to 24 cm				Adjustable from 6" to 9 ¾"/15 to 25 cm					
Seat-to-Backrest Angle	Adjustable from 95 to 115 degrees									
Maximum Tilt (Depends on Adjustments)	34 degrees				23 degrees					
Weight without Accessories	46lbs/21kg	49lbs/22kg	52lbs/24kg	62lbs/28kg	65lbs/29kg	70 lbs/32 kg	73lbs/33kg	76 lbs/34kg	86lbs/39kg	89lbs/40kg
Weight Capacity		250lbs/113kg		400 lbs	/181kg		250lbs/113kg	5	400 lbs	/181kg

Technical specifications are subject to change without notice.

Dyn-Ergo & Evolution



Applications

The S-Series Dyn-Ergo Scoot Chair and Evolution mobility chairs promote comfortable and safe independent mobility especially for individuals at risk of falls. The adjustable seat height enables efficient foot propulsion while tilt-in-space reduces sliding, promotes proper posture and helps distribute pressure. The Evolution permits both foot and/or arm propulsion and its center mounted wheelchair wheels are easily operated by users with limited arm strength and range of motion. A wide range of adjustments, sizes, and options makes configuring a Dyn-Ergo or Evolution to individual needs possible and re-assignment easier. Moreover, their friendly appearance, improved mobility and enhanced comfort make them a welcomed alternative to a geriatric chair or wheelchair.

Our products - The complete line



Decorative Metal Thera-GlideSafety Glider

Transitional Style Thera-GlideSafety Glider

Style l e® Stac Her

Utility Stacking Chair

Versatile

Versatile Adjustable Patient Room Chair



Pre-CurePressure Relief-Insert
Ready Cushion



Anti-Slip SurfacePosturo-Pedic
Positioning Aid

Contact us for more details on our range of products, or visit our website.