

Convaid Cruiser

The Convaid Cruiser is a popular 30° fixed-tilt pediatric buggy that is convenient, lightweight, and compact folding. The Cruiser is made to order to meet each individual requirements. It comes equipped with Convaid's proprietary self-tensioning seating system, Self-Tensio®, which provides increased seating tolerance and can help promote neutral pelvic positioning. It provides growth capability through its seat depth extension system.



Lightweight, Compact-Folding, and Convenient

The Convaid Cruiser is a 30° fixed-tilt buggy that improves postural control through upper body positioning.



Tension adjustable back upholstery

Adjustable back straps allow for individual adjustment to the Lumbar and Thoracic sections for the back.



One-piece height adjustable push handle

Two-piece push handle option



Self-Tensio®

Self-Tensio® is Convaid's unique self-tensioning seating system which provides pelvic positioning similar to a solid seating system in an umbrella folding rehab buggy.



Tool-less upholstery removal

No tools required to remove for cleaning

Standard

- 30° Fixed Tilt
- One Piece Height Adjustable Push Handle (CX 10/12/14/16)
- Two Piece Push Handles (CX 18)
- Tension Adjustable Back
- Tool-less Upholstery Removal
- Tool-less Adjustable Seat Depth Growth of 13mm
- Zippered Pockets in Seat & Back Upholstery (Cordura Only: CX10/12/14/16)
- 2-Point Positioning Belt
- Self-Tensio®
- Elevating Legrest
- Removable, Height Adjustable Swing-Away Footplates
- 191mm x 51mm Front, 292mm x 64mm Rear Solid Tire (CX10/12/14/16)
- Quick Release Rolko Rear Wheels (CX10/12/14/16)
- Foot Operated Wheel Lock (CX10/12/14/16)
- Toggle Wheel Locks (CX18)
- Silver Powder Coated Frame
- Made in the USA

Options

Width Sizes (mm)

- 254
- 305
- 356
- 406
- 457

Seat Depth Tubes

- 51mm seat depth tube extensions
- 127mm seat depth tube extensions
- Custom seat depth extension 76mm

Upholstery

- Cordura
- Heavy Duty Cordura

Tires

- 317mm Rear Solid Knobby Tire
- 191mm x 51mm Front, 318mm Rear Pneumatic Knobby Tire
- 203mm x 51mm Front Ribbed, 318mm Rear Solid Knobby Tire
- 203mm x 51mm Front Ribbed, 318mm Rear Pneumatic Knobby Tire

Cordura Upholstery Colors



Princess Pink



Cherry Red



Apple Green



Forest Green



Aqua Blue



Royal Blue



Sassy Purple



Chocolate Brown



Panther Black

Frame Color



Cosmic Silver

Frame Options



Scout



Planar (Not available in EU)

Accessories



Extended Headrest Cover
- w/and without windows
(Canopy)



Headrest Cover (Canopy)



Rain Cover



Mosquito Net



**Headrest Extension -
Rectangular - Textilene**

(Not available in EU)



**Headrest Extension -
Curved - Cordura**



Occi-Headwings



Padded Headwings



Seat Reducer Insert



Full Torso Support Vest



**Soft Adjustable Lateral
Support Single Flap w/
Scolio Strap**



**Soft Adjustable Lateral
Support Double Flap w/
Scolio Strap**

Accessories



H-Harness with padded covers

3-Point Positioning Belt Required



General Use Cushion*¹

(Not available in EU)



Align Cushion - Lateral Pelvis/Lateral Thigh Support*¹

(Not available in EU)



Position Cushion - Medial Thigh Support/Anti Thrust*¹

(Not available in EU)



3-Point Positioning Belt



Lateral Thigh Support - Abductor



Medial Thigh Support - Abductor



Flip-Up Armrests



Anatomic Back

Heavy Duty Reinforced Upholstery Required



Heavy-Duty Reinforced Upholstery - Cordura



Attendant Hand Brake



Upper Extremity Support Surface (Clear Tray)

[^] Crash tested Foam Tray with vinyl cover. Flip-Up Armrests required"

*Manufacturer Installation Required

¹Incontinence Cushion Liner available

Accessories



Under Seat Storage Basket - Fixed Tilt



Saddle Bags



Utility Bag




Travel Bag



Calf Panel



Angle Adjustable Footplates



Footplate Securement Strap



Castor Lock for 2" Wide Front Tires
(Standard in EU)



5-Point Harness with Padded Covers (EZ 12 only)*
Crash Tested, WC19 Compliant



Q'Straint Transit Lap Belt
(Not available in EU)



Sure-Lok Transit Lap Belt
(Not available in EU)



Transit Anchors

*Manufacturer Installation Required