

Jazzy® Carbon¹

HCPCS Code	K0899 Power Mobility Device, not coded by DME PDAC or does not meet criteria
Weight Capacity	300 lbs.
Turning Radius ³	24"
Overall Width ³	24"
Overall Length ³	37"
Maximum Speed ²	Up to 3.7 mph
Ground Clearance ³	2"
Front Casters ³	7"
Drive Wheels ³	8.5" solid
Rear Casters ³	N/A
Drivetrain	24v/250W brushless motors
Braking System	Regenerative and electromechanical
Suspension Type	Front
Standard Seat ³	17" x 17" black foam
Seat-to-Floor Range ³	20.25" with standard cushion 18.75" without cushion
Per Charge Range ^{3,5}	Up to 9.3 miles
Standard Electronics	35A Brushless Micon
Component Weights ⁶	Base: 39.6 lbs. Battery: 4 lbs.
Total Weight ⁴	43.6 lbs. with batteries
Battery Requirements	24V 12AH (288WH) lithium-ion
Battery Charger	2 amp, lithium-ion off-board
Warranty	Lifetime limited on frame; 1-year limited on electrical control/ lever system, motor/drive system, bearing sleeve and battery

- Jazzy® Power Chairs are FDA Class II Medical Devices designed to aid individuals with mobility impairments.
- Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).
- Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).
- Total weight includes base, seat and batteries.
- Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- Battery weight may vary based on manufacturer.

PRIDE
MOBILITY
LIVE YOUR BEST®

(US) 800.800.8586 • (Canada) 888.570.1113

pridemobility.com |    

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement.

The information contained herein is correct at time of publication; we reserve the right to alter specifications without prior notice.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

0004313-092922



Jazzy
carbon

TRAVEL LITE

FDA Class II
Medical Device¹

AMERICA'S FAVORITE POWER CHAIR

pridemobility.com |    

TRAVEL LITE

Padded seat back with storage pocket

Ultra-lightweight, carbon fiber frame

Airline compliant lithium battery*

Under-seat storage

Removable joystick with integrated USB charger

Compact base easily folds and unfolds for easy transport in a few easy steps

Front suspension for a smooth ride

FDA Class II Medical Device¹

Travel Lite

You know Jazzy® for its industry-leading performance and quality. Now let's talk about portability. The Jazzy® Carbon uses ultralight but exceptionally strong carbon fiber technology for its frame and a lightweight, reliable lithium-ion battery.

Travelers keep their Jazzy Carbon charged and standing by for spontaneous jet-setting adventures. They know it breezes through airport security, easily stows on the plane and gives them the mobility support they need to make the most of their travels.

Weighing in at a stunning 43.6 pounds, the Jazzy Carbon still provides top speeds up to 3.7 mph with a weight capacity of 300 pounds.

- Joystick with built in USB charger is removable for transport
- Airline compliant lithium battery*
- Front suspension for a smooth ride
- Ultra-lightweight, carbon fiber frame specifically designed for traveling
- Compact base easily folds and unfolds for easy transport in a few easy steps
- Padded seat back with storage pocket and soft, removable seat cushion for all-day comfort
- Free-wheel levers easily disengage the brake, allowing this lightweight power chair to be pushed like a manual chair
- Under-seat storage and adjustable lap belt

Available Colors



Black Carbon Fiber

Accessories

- Stealth® Simplicity cushion
- Military patches
- Additional battery

* Transport of this item on commercial aircraft is subject to federal regulations and airline policy. Please consult your airline prior to your flight to ensure your Jazzy® Carbon and its lithium ion battery are permitted.

Military Patches



Removable joystick with built in USB charger

Folds quickly and easily for transport or storage

Front suspension provides a smooth ride

Height: 28.3" | Width: 24" | Depth: 12.2"