

Convaid EZ Rider

The Convaid EZ Rider is a popular upright pediatric buggy that is convenient, lightweight, and compact folding. The EZ Rider is made to order to meet each individual requirements. It comes equipped with Convaid's proprietary self-tensioning seating system, Self-Tensio®, which provides increased seating tolerance and can help promote neutral pelvic positioning. It provides growth capability through its seat depth extension system.



Lightweight, Compact Folding, and Convenient

The Convaid EZ Rider is a 10° fixed-tilt buggy that provides upright posture to enhance functional positioning.



Tension adjustable back upholstery

Adjustable back straps allow for individual adjustment to the Lumbar and Thoracic sections for the back.



Upholstery with built-in back and seat zippers

Allow for foam inserts for additional padding and support.



Self-Tensio®

Self-tensio is Convaid's unique self-tensioning seating system which provides pelvic positioning similar to a solid seating system in an umbrella folding rehab buggy.



Upholstery with built-in seat zippers

Allow for foam inserts for additional padding and support.

Standard

- 10° Fixed Tilt
- One Piece Height Adjustable Push Handle (EZ 12/14/16)
- Two Piece Push Handles (EZ 18)
- Tension Adjustable Back
- Tool-less Upholstery Removal
- Tool-less Adjustable Seat Depth Growth of 5"
- Zippered Pockets in Seat & Back Upholstery (Cordura Only: EZ 12/14/16)
- 2-Point Positioning Belt
- Self-Tensio®
- Removable, Height Adjustable Swing-Away Footplates
- 191mm x 51mm Front, 279mm x 64mm Rear Solid Tires (EZ 12/14/16)
- 203mm Front, 279mm Rear Solid Tires (EZ 18)
- Quick Release Rolko Rear Wheels (EZ 12/14/16)
- Foot Operated Wheel Locks (EZ 12/14/16) (Toggle Wheel Locks (EZ 18))
- Silver Powder Coated Frame
- Made in the USA

Options

Width Sizes (mm)

- 305
- 356
- 406
- 457

Seat Depth Tubes

- 51mm seat depth tube extensions
- 127mm seat depth tube extensions
- Custom seat depth extension 76mm


Upholstery

- Cordura
- Heavy Duty Cordura

Tires

- 191mm x 51mm Front, 318mm Rear Pneumatic Knobby Tires
- 318mm Rear Solid Knobby Tires
- 203mm x 51mm Front Ribbed, 318mm Rear Solid or Pneumatic Knobby Tires
- 203mm x 51mm Front Solid Tires

Cordura Upholstery Colours



Princess Pink



Cherry Red



Apple Green



Forest Green



Aqua Blue




Royal Blue



Sassy Purple



Chocolate Brown



Panther Black

Frame Options



Convertible



Planar (Not available in EU)

Accessories



**Extended Headrest Cover
w/ and without windows
(Canopy)**




Headrest Cover (Canopy)



Rain Cover
Canopy Required



Mosquito Net
Canopy Required



**Headrest Extension
Rectangular - Textiene**
(Not available in EU)



**Headrest Extension
Curved - Cordura**



Occi-Headwings



Padded Headwings



Seat Reducer Insert



Full Torso Support Vest



**Soft Adjustable Lateral
Support Single Flap w/
Scolio Strap**



**Soft Adjustable Lateral
Support Double Flap w/
Scolio Strap**

Accessories



H-Harness with padded covers

3-Point Positioning Belt Required



General Use Cushion*¹

(Not available in EU)



Align Cushion - Lateral Pelvis/Lateral Thigh Support*¹

(Not available in EU)



Position Cushion - Medial Thigh Support/Anti Thrust*¹

(Not available in EU)



3-Point Positioning Belt



Lateral Thigh Support (Adductor)



Medial Thigh Support (Abductor)



Height Adjustable Flip-up Armrests



Upper Extremity Support Surface (Clear Tray)



Upper Extremity Support Surface (Foam Tray)[^]



Heavy Duty Reinforced Upholstery - Textiene

(Not available in EU)



Heavy Duty Reinforced Upholstery - Cordura

[^] Crash tested Foam Tray with vinyl cover. Flip-Up Armrests required

*Manufacturer Installation Required

¹Incontinence Cushion Liner available

Accessories



Utility Bag



Travel Bag



Travel Bag for Convertible
Wheels



Calf Panel



Individual Angle
Adjustable Footplates



Footplate Securement
Strap



Rear Anti-Tip Tubes (Pair)



Caster Locks
(Standard in EU)



5-Point Harness with
Padded Covers
(EZ 12 only)*

Crash Tested, WC19 Compliant
(Not available in EU)


*Manufacturer Installation Required



Q'Straint Transit Lap Belt
(Not available in EU)



Sure-Lok Transit Lap Belt
(Not available in EU)



Transit Anchors