

Contour™ Supreme Bath Chair

BATHING



ACCESSORIES



Rolling Shower Base



11" Extension Legs



15" Extension Legs



Head Support



Lateral Support



Positioning Belt



Designed as an alternative to our stainless steel bathing products, the Contour™ Supreme Reclining Bath Chair with the combination of ComfortMesh™ and Medical Grade PVC Base offers a comfortable seating environment with an extremely lightweight but stable foundation all without sacrificing functionality.

AVAILABLE COLORS



Beach Bubble Blue



Lionfish Lavender



Folds flat for convenient storage

BENEFITS

User Accommodation

FEATURES

- Provides parent or caregiver multiple bathing options with the simple push of a button (no tools required),
 - Five different seat to back angles can be adjusted on the fly (0°, 15°, 30°, 45°, 60° degrees) and
 - Three different leg to seat angles can be created (0°, -15°, -30°)
- Two standard positioning belts that can be used in multiple positions provide postural support.

Durable Materials

- Bacterial and fungal resistant ComfortMesh seat cover and Medical Grade PVC combine to provide a great long-term solution for daily bathing at home or in a clinical setting.

Optional Accessories

- Medical Grade PVC base available in two sizes that increase chair height by either 11" or 15"
- Additional belts and supports can be added to further enhance positioning efforts

For more information, call us today!