

EZ Lift Assist

**EZliftASSIST™**
A SHIELD INNOVATIONS PRODUCT



Only 25" wide • Holds 425lbs • 96W Electronic Motor to Raise & Lower
Locking 360 degree Casters • 12-15 day Battery Life • Rechargeable
Lithium Battery • Removable IV Drip Bag Holder • Patient Foot Rests
Optional Food Tray & Seat Cushion

Available options include
a full seat cushion, and mounted tray.



The EZ Lift Assist is a unique device designed for transporting mobility-disabled persons from bed, or wherever they may be seated, to the bathroom without having to stand or use physical strength to get there. It lifts over 400 pounds gently and easily, and rolls on floating 360-degree swivel, locking casters. The EZ Lift Assist also comes with a removable IV drip bag holder and a rechargeable, long-lasting battery. Available options include a full seat cushion, and mounted tray.

- Only 25" wide
- Holds 425lbs
- 96W electric motor
- 360-degree locking casters
- 12-15 day battery life
- Rechargeable lithium battery
- Rubber-coated foot rests
- Lifetime warranty on the frame
- 3-year warranty on the battery

Manufactured by Shield Innovations LLC