

LINCOLN EXPANDABLE BARIATRIC BED



With BUILT-IN DIGITAL SCALE SYSTEM



- **ULTIMATE PATIENT SAFETY**
- Weighs Patient in Bed
- Large Easy to Read Color Digital Display Panel

ILLUMINATED SCALE SYSTEM with Nurse Control & Patient Override Function

- Nurse Control and Pendant control all 10 positions
- Lockout feature on nurse control overrides patient hand pendant
- Large illuminated color screen displays patient weight
- Stores weight history and calculates difference in weight
- LB./KG selectable
- CPR and emergency stop function
- Comes pre-calibrated
- Powered by AC or backup 12V rechargeable battery included

MODEL # **LX-BARI-S**



DIGITAL SCALE / NURSE CONTROL FEATURES & FUNCTIONS

Digital Scale System - weighs patient in the bed

2 white LED's under the nurse control illuminate the foot end at night

Bed alarm function - sounds when patient is off the bed

Battery power indicator on display

LINCOLN EXPANDABLE BARIATRIC BED

WITH SPLIT FRAME DESIGN

Ideal Rental Solution - folds compact for easy transport and storage



MODEL # **LX-BARI-SF**

| Model # | Width | Length | Low Range | High Range | Weight Cap |
|----------------------------|---------------|----------|-----------|------------|------------|
| LX-BARI | 36"- 42"- 48" | 80"- 88" | 14" | 30" | 750 |
| LX-BARI-9 (Ultra Low) | 36"- 42"- 48" | 80"- 88" | 9" | 25" | 750 |
| LX-BARI-S (Built-In Scale) | 36"- 42"- 48" | 80"- 88" | 14" | 30" | 750 |
| LX-BARI-SF | 36"- 42"- 48" | 80"- 88" | 14" | 30" | 750 |