

# POWERED PATIENT LIFT BARIATRIC WITH DUAL MOTORS



## FREE:SP:RIT

Model: FS660 Bariatric Patient Lift



Bariatric

## PROVIDES A SAFE ALTERNATIVE TO MANUAL LIFTING FOR BOTH THE CAREGIVER AND PATIENT

### SECURE AND STURDY



#### ● DUAL MOTORS

Raises/Lowers the patient from floor to 70" high

#### ● POWER BASE

Expands from 28" - 45" ultra wide

**660 Lbs.**  
Weight  
Capacity



- Dual motors power the lift range and legs to open and close
- Ultra wide leg opening easily wraps around bariatric wheelchairs and equipment
- 50" extra long leg length for increased stability
- Oversized EZ Grip foam handle for increased control when moving the lift
- Emergency stop button locks the lift in place
- CPR function lowers and raises the lift when the battery is discharged
- Total pinch point protection
- Padded swivel bar protects residents head



- Quick release battery eliminates downtime
- Battery can be replaced immediately with separate charging base
- Battery can also be charged on board
- LED's show battery charge level
- Charger case stores the charger plug on the lift

**WEIGHT CAPACITY:** 660 Lbs.

**LIFTING RANGE:** 28" - 70" high

**BASE WIDTH:** Open- 45" / Closed- 28"

**BED CLEARANCE:** At least 4.5"

**CASTER SIZE:** Rear- 4" locking / Front- 3" swivel

# POWERED PATIENT LIFT BARIATRIC WITH DUAL MOTORS



## FREEspirit

Model: FS660 Bariatric Patient Lift



Bariatric



**Power Base**  
Adjustable Width from 28" - 45"



**3 Oversized locking hooks**  
on each side  
Safely secures sling to lift



**Dual Function Controls**  
- by remote handset  
or on board

**CPR EMERGENCY  
STOP & LOWERING FUNCTION**



**Bariatric Digital Scale**  
(Optional)

Large easy to read  
display with bright  
blue backlight



**Charger Plug Storage  
Case - Included**



**Front and Rear  
Locking Casters**  
Secures Lift in Place

MODEL # **FS660-T**