

ROLLING SHOWER WHEELCHAIR



ROLLING SHOWER CHAIR WITH DETACHABLE COMMUNE

SHOWER PROOF WITH SEALED GASKETS



- Castors: Front: 2 x 5" / Rear: 2 x 24"
- Seat width: 18"
- Seat depth: 16"
- Seat height: 20" - 22"
- Back height: 17"
- Height of arm: 9"
- Total width: 27"
- Total height: 38" - 40"
- Total length incl. leg rests: 40"
- Weight capacity: 265 lbs.
- Weight: 40 lbs.



5" Swivel wheels with foot pedal locking castors



Commune can be easily removed without lifting patient or raising the seat



Padded seat with a removable section for easy cleaning



24" wheels with hand-brake allow patients to maneuver the chair without an aid



Arm drops for easy transfer



Detachable and height adjustable footrests

MODEL # **SWC1816**