

# BARIATRIC COMMODE



MODEL # **CB1000**

- **WIDTH**  
25.75" between arms /  
26.25" actual width
- **SEAT DIMENSIONS**  
23.5" W x 18.5" D
- **HEIGHT**  
17.5" - 22"
- **WIEGHT**  
21 LBS.
- **WEIGHT CAPACITY**  
1,000 LBS.



**Bariatric**

- Easy to release arm mechanism allows for safe lateral patient transfers to and from commode
- Durable and heavy-duty steel tubing is safe for larger patients
- Large durable snap on seat

- Padded armrests provide comfort and stability
- Pail easily slides in from front of the commode
- Comes complete with 12 qt. commode bucket and cover that removes from the front
- Available in heavy and ultra heavy weight capacities

- **WIDTH**  
24" between arms / 32" actual width
- **SEAT DIMENSIONS**  
14" W x 17" D
- **HEIGHT**  
16" - 22"
- **WIEGHT**  
19 LBS.
- **WEIGHT CAPACITY**  
650 LBS.



MODEL # **CB650**