

MATTRESS Standard Size



ALTERNATING PRESSURE COMBINED WITH LOW AIR LOSS

Effective in Preventing
Stage IV Wounds

POLY/PU QUILTED COVER
FOR ULTIMATE IN PATIENT COMFORT



Whisper Quiet
Digital Pump



MODEL # CZ36



Therapeutic Support Mattress Replacement System
Designed to prevent, treat and manage pressure ulcers



Alternating Pressure with Low Air Loss
Cell-on-cell Design Prevents Bottoming Out

MATTRESS Standard Size



PUMP



Alternating Pressure
and Static Function




Comfort Level
Adjustable



Low Pressure
Indicator



Panel Lock
Function



- Auto Firm Inflation
- Low Pressure Alert
- Power Failure Alert
- Quick-change Air Filter

USER FREINDLY DIGITAL CONTROL PANEL

PRESSURE RANGE - 10 LEVELS
Adjustable by patients weight
or pressure range

CYCLE TIME VARIABLE:
10, 15, 20, 25 minutes

AIR OUTPUT:
8 LPM

PUMP DIMENSIONS:
11.4" x 12" x 5.9"

MATTRESS



20 Easy Replaceable
snap out Air Cells



Cell-on-Cell
Design

- Pillow Function - Top 3 cells remain static for patient comfort
- Patented CPR Quick Inflate/Deflate System
- 8 Ventilated Low Air Loss Cells



EZ SNAP-IN-SNAP-OUT
REPLACABLE CELLS

CELL MATERIAL:
Nylon/TPU

CELL NUMBER:
20 Cells

COVER MATERIAL:
Fluid Resistant Nylon & Quilted

BASE MATERIAL:
Woven with non-slip material

MATTRESS DIMENSIONS:
80" L X 36" W X 8" H

WEIGHT CAPACITY:
450 Lbs.

Also Available with a
Built-In Air Perimeter

Prevents Falls and for
Ultimate Patient Safety



Model #	Description	Dimensions	Alternating Pressure	Low Air Loss	Cell-On-Cell	Weight Capacity
CZ36	Standard Length	80" X 36" X 8"	Yes	Yes	Yes	450 Lbs.
CZ3684	Extra Long	84" X 36" X 8"	Yes	Yes	Yes	450 Lbs.
CZ36-BIP	Built-In Perimeter	80" X 36" X 8"	Yes	Yes	Yes	450 Lbs.