



### High Specification Self Adjusting Mattress System

- Choose the most appropriate therapy mode to suit individual needs of patients.



### Optional: Alternating Pressure Pump

- Non-powered mattress replacement system featuring self adjusting immersion therapy.
- Powered optional alternating pressure control unit for high-risk of pressure injury.
- Superior four(4)way stretch micro climate cover
  - Low friction and shear
  - immersion fabric is fluid-proof and vapor permeable.




## AirPlus Self Adjusting Immersion Mattress

### Deep air chamber, self adjusting immersion therapy.

Patented self-adjusting technology allows internal cells to automatically adjust to each individual body type and weight.

Head and back section engineered to provide torso and back support.

**Lower extremity section** molds to the lower leg to accommodate limb contours and off-load heel/ankle skeletal structures; areas prone to pressure injury development.



### SAFETY • IMMERSION • MOBILITY

• Edge of mattress support bolster can assist with sitting and standing balance upon bed exit.

-Promotes roll stability for fall prevention.

#### Sizes Part Number

36" x 80"	- MPAM3680SAT
36" x 84"	- MPAM3684SAT
39" x 80"	- MPAM3980SAT
39" x 84"	- MPAM3984SAT
42" x 80"	- MPAM4280SAT
42" x 84"	- MPAM4284SAT
48" x 80"	- MPAM4880SAT