

Battery-Powered Patient Lift

HCPCS: E0635



- Heavy-gauge construction
- Jumbo actuator pump with emergency button can switch to manual mode to lower patient down safely
- Fast, easy charge from AC outlet
- To charge, simply plug in unit Battery does not need to be removed to be charged Lift will not operate while plugged in For home-care use
- Audible warning when battery power falls below 50%

PRODUCTS	Description
13240	Patient Lift - Electric, 1/cs
SPECIFICATIONS	Description
1 Year Electronics	Battery, Hand control, Actuator
Base Clearance	2"
Base Height	4.75"
Base Width	41" (open) 24" (closed)
Boom Height	27.5" - 76"
Carton Shipping Weight	122 lbs.
Casters	5" Rear
Motor	Linak Motor
Warranty	Limited Lifetime
Weight	108 lbs.
Weight Capacity	450 lbs
PARTS	Description
13240BATT	BattPackw/Batteries,13240,13244,13246,1/, EA
13255	External Charger w/out PowerCord,1/ea, EA
13260	Battery&ChargerBracket,Elec Patient,1/ea, EA
13255PC	Charger Power Cord,13246,1ea, EA
13240AP	Spreader Locking Plate,13240,1ea, EA
13240A	Actuator for 13240,1/ea, EA
13240SBH	Handgrip,SpreaderBar,13240,1/ea, EA
13240CA	Caster w/o Lock,13240,1/ea, EA
13240FP	Square Frame Plug,13240,1/ea, EA
13240BC	Battery Cell only for 13240,1/ea, EA
13240HC	Hand Control for 13240,13244,13246,1/ea, EA
13240PC	Power Cord for 13240,13244,1/ea, EA

PARTS	Description
13240ARB-N	ActuatRodBoltWithMetalSpacer, 1/pkg, PKG
13240-CC-SL	Cradle Clips, Spring-Loaded,13240, 1/pr, PR
CPCFUSE	Fuse,10A,250V,CPC,13240,1/ea, EA
13240CAB	Cradle Attchmnt Bolt,13240,1/ea, EA
13240MU	Mast, Upper,13240,1/ea, EA
13058PB	Cradle Pivot Bolt,13010SV,1ea, EA
13240-CC	Cradle Clips for 13240, 1/pr, PKG
13240PH	Push Handle for 13240, 1set, EA