

Pride®

Lifts & Ramps

PRIDE

TOUGH

The Milford Person Lift is a highly adaptable and versatile wheelchair to vehicle transfer solution renowned to be quick and easy to operate for a wide range of people. With a lifting capacity of 330 lbs., the Milford Person Lift can be used anywhere from transferring to your vehicle, in the home and office, or while traveling.

FEATURES

- 330 lbs. weight capacity
- Works in a wide variety of vehicles
- Can be fitted to the driver's side, passenger's side, and front, middle or back seats
- Sling sizes from extra small to extra large with amputee slings also available
- Lifts quickly and easily with the touch of a button
- Stress-free application for the patient and caregiver
- Emergency switch and optional manual operation

Milford Person Lift

BASE WEIGHT ¹	18 lbs.
WEIGHT CAPACITY	330 lbs. max.
MAXIMUM TRAVEL DISTANCE FROM THE GROUND	22.5"
MOTOR	12 Volt DC
WARRANTY	1-Year Transferable
HITCH CLASS	N/A
OPTIONS	Sling Size & Type, Wall Bracket, Smartbase
VEHICLE PLATFORM	Full-size van, Minivan, Sedan, SUV and Crossover
LIFTS CURRENT PRIDE® PRODUCTS	N/A

1. Does not include any accessories or hardware.

Truly adjustable transfer solutions, quick & easy to operate



Pride®
Mobility Products Corp.
— LIVE YOUR BEST! —

prideliftsandramps.com • pridemobility.com
(800) 800-8586 (U.S.) • (888) 570-1113 (Canada)

Milford SmartBase

Smartbase option available for mobile application

OPTION

MILFORD-FY-12-02-15

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.