

PreserveTech Secure Lock Raised Toilet Seat



HCPCS: E0244



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F



FIGURE G

- Three locking points of security located on the front, left and right side providing the highest level of stability and security. Most raised toilet seats only offer one or two locking points of security.
- Adding 5 inches of height to your toilet making it easier and safer to sit or stand.
- Padded armrests provide superior support and create a comfortable, non-slip surface for a safe and secure grip when sitting or standing up.
- Easily adjust armrest height to your needs.
- Contoured to provide a safe and comfortable elevation.
- Made from HDPE that is both lightweight and durable to support up to 350 lbs.
- Fits comfortably on most standard and elongated toilets.
- No Tools Required! Simply turn the adjustable knobs to tighten the locking plate and ensure a secure fit.
- Made using recycled materials.
- Drive PreserveTech products are manufactured with an antimicrobial* agent that helps to minimize the growth of stain and odor-causing bacteria and mold on treated surfaces.
- *Antimicrobial properties are built-in to protect the products. These products do not protect users or others against bacteria, viruses, germs or other disease organisms.
-

PRODUCTS	Description
RTL12C003-WH	PreserveTech Secure Lock RTS, 1EA, 1/CV
SPECIFICATIONS	Description
Product Type	Toilet Seat
Seat Depth	20.3"
Seat Width	17.5"
Waranty	Limited Lifetime
Weight Capacity	350 lbs
Width Between Arms	17.5"
Material Composition	Seat: HDPE Armrests: Aluminum