

2-in-1 Locking Raised Toilet Seat with Tool-free Removable Arms

HCPCS: E0244



FIGURE A



FIGURE B



FIGURE C

- Designed for individuals who have difficulty sitting down or getting up from the toilet
- Adds 5" to height of toilet seat
- Tool-free arms can be removed or added (Figure A)
- Arm sleeves and receivers are made of metal
- Heavy-duty molded plastic construction
- Locking device with larger, heavy-duty "worm screw" and locking plate for a safe, secure fit
- Lightweight and portable
- Easily attaches to toilet bowl, tool-free
- Easy-to-clean
- Fits most toilets
- Model RTL12027RA is retail packaged
- Made using recycled materials

PRODUCTS	Description
12027RA-4BULK	Raised Toilet Seat w/RemArms,Bulk,4/cs
RTL12027RA	Elev.Toilet Seat w/ RemArms,RTL,1 EA
SPECIFICATIONS	Description
Carton Shipping Weight	8.5 23 lbs.
Seat Dimensions	17" (W) x 16.5" (D) x 5" (H)
Warranty	Limited Lifetime
Weight	5 lbs.
Weight Capacity	300 lbs
Width Between Arms	18"