

Bariatric Folding Commode

HCPCS: E0168



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E

- Tool-free assembly in 3 easy steps (Figure A - C)
- Easily opens and folds for storage and transportation
- Durable plastic 16 ½" deep snap-on seat with lid
- Plastic armrests for added comfort
- Removable tool-free back
- Comes complete with 12 qt commode bucket with carry handle, cover and splash shield
- Weight Capacity 650 lbs
- Made using recycled materials

PRODUCTS

Description

11117N-1	FoldingBariatric Commode (DeepSeat),1/cs
----------	--

SPECIFICATIONS

Description

Carton Shipping Weight	25/50 lbs.
Height	15.5" - 22"
Seat Dimensions	13.75" (W) x 16.5" (D)
Warranty	Limited Lifetime
Weight	19 lbs.
Weight Capacity	650 lbs
Width	24" between arms 31.75" Outside