

Competitive Edge Line 3-in-1 Folding Commode

HCPCS: E0163



- Easily opens and folds
- Folds flat for convenient storage and transportation
- Tool-free set-up
- Durable plastic snap on seat
- Grey powder coated steel welded construction increases strength and durability
- Comes complete with 7.5 qt commode bucket
- Splash Guard included
- Limited Lifetime Warranty
- Made using recycled materials

PRODUCTS

Description

11148CE-4	Commode,Folding ComBid,4/cs
-----------	-----------------------------

SPECIFICATIONS

Description

Carton Shipping Weight	61 lbs.
Seat Dimensions	13.5" (W) x 15" (D)
Seat Height	16.6" - 22.5"
Warranty	Limited Lifetime
Weight	15 lbs.
Weight Capacity	350 lbs
Width	18" Between Arms 22.5" Outside