

Combination Padded Transfer Bench/Commode

HCPCS: E0247



FIGURE A



FIGURE B



FIGURE C



FIGURE D

- Combines transfer bench and commode into one
- Tool-free assembly of back, legs and arm (Figure A - C)
- Comfortable cushioned seat and backrest
- Aluminum A-frame construction is lightweight and provides additional stability
- Height adjusts in ½" increments with unique "Dual Column" extension legs
- No exposed hardware to injure patient
- Extra-large suction cups provide added safety
- New pinch-free cover allows for push pins to be depressed without pinching finger (Figure D)

PRODUCTS

Description

12005KDC-1	TranfBnchw/CommodePaddedSeatKD,1/cs
------------	-------------------------------------

SPECIFICATIONS

Description

Carton Shipping Weight	15 lbs.
Outside Legs	17.25" (W) x 28.5" (D)
Seat Dimensions	24" (W) x 16" (D)
Seat Height	18.5" - 22.5"
Warranty	Limited Lifetime
Weight	10 lbs.
Weight Capacity	400 lbs

ACCESSORIES

Description

12011SL	SuctionCups,Lg.12011,1pr/bx, BX
12011SS	SuctionCups,Sm.12011,1pr/bx, BX