



A Division of United Care B.V.




Monarch Portable Ceiling Hoist

The Monarch Portable Ceiling Hoist is one of the lightest on the market weighing only 4.8kg, whilst lifting the industry standard of 200kg. Ideal for occasional or multiple room use, the product is robust and charges with an AC Charger. The unique spreader bar includes our industry leading Secure Sling Mechanism designed to ensure loop attachments cannot be unintentionally removed.

Perfect for rental properties, multi-room, temporary or short-term use, the Monarch Portable Ceiling Hoist can be used as a flexible alternative to fixed overhead hoists.

Unique Features

- Ultra-Lightweight, Weighing only 4.8kg
- Increased Battery Capacity
- Carry handle
- Express Recharge
- Unique Spreader Bar with Secure Sling Mechanism
- 130kg/200kg safe working load
- Ideal for transporting with the user
- Ideal for rental properties or hotels for nonpermanent use



Easy-to-use Carabiner clip



4X more powerful



50% more lift cycles



Available with telescopic reacher



Easy Transport Handle



Unique Secure Sling Mechanism

Product Specifications Monarch Portable Hoist

Specifications

Model	Monarch Portable Ceiling Hoist
Function	Patient Lifting, Room-to-Room Transfer, Temporary Use, Rental Properties
Total Weight	4.8 kg
Weight Capacity	130 kg (Standard) / 200 kg (Optional)
Lifts per Charge	30 to 40 cycles
Speed with 200 kg	4.5 cm/s
Travel	2300 mm
Soft Start / Stop	Yes
Emergency Stop / Lower	Yes
Options	Standard Portable Trolley (OpeMed track compatible), Bespoke Portable Trolley
Track Options	(compatible with leading track manufacturers), Reacher. H system, Turntable, Gate, Curves, Room-to-Room, Freestanding Gantry : 2m, 4.5m, 6m, 8m
Track Compatibility	6m, 8m
Application	Compatible with all OpeMed track and other leading track manufacturers Nursing Home, Hospital, Domestic, School, Leisure Centre, Physiotherapy

Dimensions

