



Zippie Voyage

February 2023

EI218 (With Seating - E1232)

Date: _____ Dealer Acct #: _____ Dealer: _____ Dealer Contact: _____ Dealer Address: _____ Dealer City: _____ ST: _____ ZIP: _____ Dealer Phone: () _____ Fax: () _____ Confirm Via: ☐ Fax ☐ Email Mark For: _____	SUBMITTING FOR: ☐ Quote ☐ Order PO#: _____ Ship To: _____ Attention: _____ Address: _____ Address: _____ Ship To City: _____ ST: _____ ZIP: _____ Ship To Phone: () _____ Fax: () _____
--	---

The **HCPCS CODES** provided are based on code verification through PDAC or our interpretation of the code definitions.
 The use of **HCPCS CODES** does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price
1. FRAME STYLE									

Designed with busy families in mind. The stroller frame quickly and easily folds after removing the seat in one simple step for compact storage. Multiple seating configurations allows you to personalize your Zippie Voyage to perfectly fit your needs. Each unit includes:

- | | | |
|---|--|--|
| <ul style="list-style-type: none"> Foot Operated Wheel Lock 12" Rear Tires & 8" Front Casters Front Caster Suspension & Swivel Locks 75 lbs Weight Capacity | <ul style="list-style-type: none"> Quick Release Wheels Angle Adjustable Stroller Handle Adjustable Seating Shell Canopy w/ breathable window/mesh pockets | <ul style="list-style-type: none"> 35° Tilt Fwd & Rwd facing seating capability Under carriage storage (may be limited when some options chosen) |
|---|--|--|

Voyage Frame Seating Options			
☐	005BS7	E1232	Voyage for Advanced Seating with Tilt - Includes Frame and Seating
☐	005BS9	E1232	Voyage for Moderate Seating with Tilt - Includes Frame and Seating

Seating Only Options (Omit Frame)			
☐	005BS15		Advanced Seating with Tilt - Omits Frame, Seat-Only
☐	005BS16		Moderate Seating with Tilt - Omits Frame, Seat-Only

Booster Base Seating Options			
☐	005BS10	E1399	Voyage Booster Base with Tilt and Advanced Seating
☐	005BS11	E1399	Voyage Booster Base with Tilt and Moderate Seating
☐	005BS12	E1399	JCM Booster Base with Tilt and Advanced Seating
☐	005BS13	E1399	JCM Booster Base with Tilt and Moderate Seating

Transit			
☐	005KJ11	K0108	Transit Option
			\$328



2. MODERATE SEATING

Choose Moderate Seating options, prices reflect ordering base which includes seating.

Moderate Seating Options

☐	005KM01		Moderate Seating	NC
☐	005KM02	E1014	Recline Feature (required)	\$391

Grab Bar (Optional)

☐	005KR01	K0108	Grab Bar	\$138
--------------------------	---------	-------	----------	-------

Choose Upholstery Color (includes inner covers and crossover canopy)

☐	005KN01		Blue	NC
☐	005KN02		Red	NC
☐	005KN03		Beige	NC
☐	005KN04		Grey	NC

2A. POSITIONING & SUPPORTS - MODERATE SEAT ONLY

Head Supports (Optional)

☐	005VJ01	E0955	Head Supports - Small (Recommended)	\$168
☐	005VJ02	E0955	Head Supports - Large	\$168

Positioning - Medial Support (Optional)

☐	005KQ01	E0957	Medial Support - Tall	\$151
☐	005KQ02	E0957	Medial Support - Contour	\$157

Positioning - Laterals (Optional)

☐	005KO01	E0956	Laterals Small - (R & L)	\$168
☐	005KO02	E0956	Laterals Large - (R & L)	\$168

Positioning - Hip Guides (Optional)

☐	005KP01	E0956	Hip Guides - (R & L)	\$196
☐	005KP03	E0956	Infant Insert (Hip Guides Req'd.)	\$141

JAY Anterior Pelvic Belts - E0978

Standard 2 Point Belt - Push Button 1"

☐	005ZZ10	E0978	No Pad (APB5110)	\$46
☐	005ZZ11	E0978	4.5" Pad (APB5115)	\$102
☐	005ZZ12	E0978	5.5" Pad (APB5116)	\$102

Standard 2 Point Belt - Push Button 1 1/2"

☐	005ZZ14	E0978	No Pad (APB7210)	\$46
☐	005ZZ15	E0978	7.5" Pad (APB7212)	\$97

V-Contour Pelvic 4 Point Belt - Push Button Varies

☐	005ZZ24	E0978	SM Pad - 1" Push Button (APB3716)	\$189
☐	005ZZ23	E0978	Med Pad - 1-1/2" Push Button (APB3727)	\$189

JAY Anterior Trunk Supports - E0960

Classic Structure not available with Early Intervention Size

Anterior Trunk Support - Style

☐	005JY01	E0960	Classic Dynamic (ATS8100 Series)	\$170
☐	005JY02	E0960	Classic Structure (ATS8300 Series)	\$170
☐	005JY06	E0960	Center Opening Dynamic (ATS8800 Series)	\$170
☐	005JY05	E0960	Center Opening Structure (ATS9400 Series)	\$170

Anterior Trunk Support - Size

☐	005JZ05		Early Intervention	NC
☐	005JZ06		Xtra Small	NC
☐	005JZ01		Small	NC

2B. FOOTREST - MODERATE SEAT ONLY

☐	005F329	K0040	Angle Adj. Footrest (Required)	\$173
--------------------------	---------	-------	--------------------------------	-------

					
Grab Bar	Head Supports	Medial Support	Laterals	Hip Guides	Infant Insert

3. ADVANCED SEATING

Choose Advanced Seating options, prices reflect ordering base which includes seating.

Planar Seat & Contoured Seat Cushion Include:

- Contoured Medium Density Foam
- 6" to 13" Hip Width
- Built in Depth Growth
- Removable Neoprene Cover
- Optional 3" Lateral Hip Blocks (Pr)
- Optional Slight Medial Thigh Contour
- Optional Fixed Ischial Ledge



Planar & Contoured Seat Options:


☐	005CU28		Planar Seat without Hip Blocks	NC
☐	005CU41	E2292	Planar Seat w/ Hip Block - Pair	\$296
☐	005CU42	E2294	Contoured Wedge Seat	\$296
☐	005CU43	E2294	Contoured Anti-Thrust (Ledge) Seat	\$296
☐	005CU2		OMIT Cushion	NC

Choose Depth (Based on 1" T-Back)

☐	005SD78		Short 5" - 8" Depth	NC
☐	005SD79		Medium 7" - 10" Depth	NC
☐	005SD80		Long 9" - 13" Depth	NC

ADJUSTABLE CONTOUR SEAT includes:

- Choice of foam density
- Adjustable Pre-Ischial Ledge
- Adjustable 1" High Foam Wedges for:
Lateral Hip and Thigh Support
Medial Knee and Thigh Support
- 3" Seat Depth Growth Kit
- Removable Neoprene Cover



Adjustable Contour Seat Options:

☐	005ST65	E2294	Adjustable Contour Seat	\$357
☐	005CU2		OMIT Cushion	NC

Choose Foam

☐	005FM1		Soft Foam	NC
☐	005FM2		Medium/Soft Foam	NC
☐	005FM3		Medium Foam	NC
☐	005FM8		OMIT Foam	NC

Choose Depth (Based on 1" T-Back)

☐	005SD78		Short 5" - 8"	NC
☐	005SD79		Medium 7" - 10"	NC
☐	005SD80		Long 9" - 13" (4" growth)	NC

3A. BACKREST - ADVANCED SEATING ONLY

BACKREST SHELL

☐	005BT111		Short 7" to 13" Height	NC											
Desired Backrest Height Setting (<i>Bold indicates default, when not picked</i>)															
☐	7"	☐	8"	☐	☒	10"	☐	11"	☐	12"	☐	13"			
☐	005BT112		Tall 10" to 17" Height	NC											
Desired Backrest Height Setting (<i>Bold indicates default, when not picked</i>)															
☐	10"	☐	11"	☐	12"	☒	13"	☐	14"	☐	15"	☐	16"	☐	17"

BACKREST OPTION

☐	005KJ9	E1014	Recline Feature (90° - 135° range)	\$376
--------------------------	--------	-------	-------------------------------------	-------

BACKREST CUSHION OPTIONS

2.5" "T" Back is used for seat depths as small as 3.5". Backrest has removable foam for additional depth growth without adjusting seat pan.

1" Thick T-Back Cushion Options:

☐	005BT35	E2291	Short 8" to 10"	NC
☐	005BT28	E2291	Medium 10" to 12"	NC
☐	005BT29	E2291	Tall 12" to 14"	NC
☐	005BT56	E2291	Extra Tall 14" to 16" (N/A w/ Short Backrest Shell)	NC
☐	005BT133		OMIT Back Cushion	NC

2.5" Thick T-Back Cushion Options:

☐	005BT115	E2291	Short 8" to 10"	\$75
☐	005BT116	E2291	Medium 10" to 12"	\$75
☐	005BT117	E2291	Tall 12" to 14"	\$75
☐	005BT118	E2291	Extra Tall 14" to 16" (N/A w/ Short Backrest Shell)	\$75
☐	005BT133		OMIT Back Cushion	NC

3B. HEAD SUPPORT - ADVANCED SEATING ONLY

Hardware allows the headrest to be flush mounted to backrest when above shell. When choosing a Headrest option, Hardware Bracket must be chosen, Hardware also available individually. Headrest Ht. Adj. Range from seat pan based on selected Backshell Height. (i.e. Short Ht Ranges 9" - 19.5" & Tall Ht Ranges 13" - 23-1/2")

☐	005HE40	E0955	Whitmyer 8" Plush	\$242
☐	005HE41	E0955	Whitmyer 10" Plush	\$242
☐	005HE42	E0955	Whitmyer SM Contour Cradle	\$242
☐	005HE43	E0955	Whitmyer Med Contour Cradle	\$242
☐	005HE29	E0955	Whitmyer Sm Adjust-A-Plush	\$289
☐	005HQ6	K0108	Whitmyer DFS (Forehead Strap) Squared for Single pad, must order Whitmyer headrest	\$128
☐	005HE31	E1028	Hardware Bracket Only (can be chosen w/o headrest)	\$113

3C. POSITIONING & SUPPORTS- ADVANCED SEATING ONLY

For the supports, Voyage will be set up at factory. Insert dimensions on last page spec. drawing. (Not required, will default to std. settings)

LATERAL & ANTERIOR TRUNK SUPPORTS

Lateral brackets height adjustment range based on Backrest Shell Height (For Example: Short Laterals adj. 3" - 12" and Tall Laterals adj. 3" - 17")

Swing Away Laterals (Adjustable 1"-8" widths)

☐	005KE41		Swing Away Bracket Adj. - Right	\$242
☐	005KE42		Swing Away Bracket Adj. - Left	\$242
☐	005KD15	E0956	Large Straight Pad 3" x 6" - Right	\$107
☐	005KD19	E0956	Large Straight Pad 3" x 6" - Left	\$107
☐	005KD16	E0956	Large Curved Pad 3" x 6" - Right	\$107
☐	005KD20	E0956	Large Curved Pad 3" x 6" - Left	\$107
☐	005KD13	E0956	Small Straight Pad 2" x 6" - Right	\$107
☐	005KD17	E0956	Small Straight Pad 2" x 6" - Left	\$107
☐	005KD14	E0956	Small Curved Pad 2" x 6" - Right	\$107
☐	005KD18	E0956	Small Curved Pad 2" x 6" - Left	\$107

JAY Anterior Trunk Support

☐	005A827	E0960	JAY Contour Dynamic EI Pad	\$163
☐	005A828	E0960	JAY Contour Dynamic XSmall Pad	\$163
☐	005A829	E0960	JAY Contour Dynamic Small Pad	\$163

Anterior Trunk H Strap Support

☐	005SS2	E0960	XX Small H Strap 5 Point	\$177
--------------------------	--------	-------	--------------------------	-------

Fixed Laterals (Adjustable 1"-10" widths)

☐	005KE24		Fixed Angle Adjustable - Right	NC
☐	005KE25		Fixed Angle Adjustable - Left	NC
☐	005KD15	E0956	Large Straight Pad 3" x 6" - Right	\$107
☐	005KD19	E0956	Large Straight Pad 3" x 6" - Left	\$107
☐	005KD13	E0956	Small Straight Pad 2" x 6" - Right	\$107
☐	005KD17	E0956	Small Straight Pad 2" x 6" - Left	\$107

Anterior Trunk Butterfly Support

Dimensions are seat to shoulder height.

☐	005KI1	E0960	XX Small Butterfly Pad (< 10")	\$118
☐	005KI2	E0960	X Small Butterfly Pad (10" - 12")	\$118
☐	005KI3	E0960	Small Butterfly Pad (12" - 14")	\$118
☐	005KI4	E0960	Medium Butterfly Pad (14" - 16")	\$118
☐	005KI5	E0960	Large Butterfly Pad (>16")	\$118
☐	005SS1	E0960	Shoulder Strap Pads (Only Available with Butterfly Support)	\$27

LOWER BODY SUPPORTS**Lateral Hip Support (Adjusts from 7-1/2" to 11" width)**

Only available with Planar seat without Hip Blocks and Adj. Contour seat.

☐	005KC6	E0956	SM (4"x 4", 3/8" thick) Right	\$107
☐	005KC4	E0956	SM (4"x 4", 3/8" thick) Left	\$107
☐	005KC5	E0956	LG (5-3/8"x 4-3/4", 1/2" thick) Right	\$107
☐	005KC3	E0956	LG (5-3/8"x 4-3/4", 1/2" thick) Left	\$107

Lateral Thigh Support (Adjusts from 3-1/2" to 11" width)

Only available with Planar seat without Hip Blocks and Adj. Contour seat.

☐	005KF6	E0953	SM (4"x 4", 3/8" thick) Right	\$107
☐	005KF4	E0953	SM (4"x 4", 3/8" thick) Left	\$107
☐	005KF5	E0953	LG (5-3/8"x 4-3/4", 1/2" thick) Right	\$107
☐	005KF3	E0953	LG (5-3/8"x 4-3/4", 1/2" thick) Left	\$107

Pelvic Positioning Belt

Must use when ordering Anterior Trunk Supports.

☐	005A14	E0978	1" Padded	\$48
☐	005A15	E0978	1-1/2" Padded	\$48
☐	005A263	E0978	1" Padded Auto Style (Dual Pull)	\$59
☐	005A256	E0978	1-1/2" Padded Auto Style	\$59
☐	005A830	E0978	JAY Medial 4-Point Belt XS Pad	\$102
☐	005A831	E0978	JAY Medial 4-Point Belt Sml Pad	\$102

Medial Support

☐	005KH12	E0957	Small (3"Hx 2-1/2"Wx 4-1/2"L)	\$135
☐	005KH13	E0957	Large (3-3/4"Hx 4"Wx 6"L)	\$135
☐	005HE77	E1028	Hardware Medial Support Brkt	\$242

3D. UPHOLSTERY COLORS - ADVANCED SEATING ONLY**Color Options - Solid**

☐	005CL9		Purple	NC
☐	005CL20		Pink Sorbet	NC
☐	005CL11		Green	NC
☐	005CL12		Navy Blue	NC
☐	005CL13		Orange	NC
☐	005CL17		Turquoise	NC
☐	005CL16		Black	NC
☐	005CL18		Burgundy	NC

Color Options - Print (Std. Trim listed behind each print)

☐	005CL19		Leopard (Std. w/ Black Trim)	NC
☐	005CL14		Butterflies (Std. w/ Purple Trim)	NC
☐	005CL15		Red Race Flags (Std. w/ Red Trim)	NC

Personalized Option (Only available with Prints)

☐	005C001		Personalized Color Selection: Choose & Enter from Solid Color Options:	\$75
--------------------------	---------	--	---	------

Canopy Pattern

☐	005UF01		Cross Over	NC
--------------------------	---------	--	------------	----

3E. FOOTREST - ADVANCED SEATING ONLY

Footrest adjustment range is 1-1/2" to 12". Choose desired preset range, if not selected Long range is Default. For the Short Range, the footplate will be inverted.

☐	005F312	K0040	Angle Adj. Footrest (Required)	\$166
☐	005F313	K0037	Preset in short range 1-1/2"-4"	NC
☐	005F314		Preset in long range 4-1/2" - 12"	NC

4. ACCESSORIES & OPTIONS

Vent & Battery Trays

☐	005LF1	E1029	Dovetail Vertical Frame Mounted Hardware for LTV ventilators <i>(limit 15lbs & includes 1 Utility hook)</i>	\$403
☐	005LF11	E1029	Dovetail Vertical Frame Mounted Hardware for Trilogy ventilators <i>(limit 15lbs & includes 1 Utility hook) (Only includes female dovetail mount)</i>	\$403
☐	005LF2	E1029	Battery and Vent/Utility Underframe Tray <i>(limit 50lbs, ideal for U1 or Grp 22 Batteries & includes 1 Utility hook)</i>	\$440
☐	005LF3	E1029	Front Mount Utility platform <i>(limit 50lbs)</i>	\$403
☐	005LF4	E1029	Rear Add on Accessory platform <i>(limit 45lbs or one Grp 24 battery & includes 1 Utility hook. Must be removed or folded before folding base)</i>	\$188

Floor Sitter/ Feeder Base

☐	005A083	E1399	Folding Booster Base(Ht. Adj.)	\$1,498
☐	005A082	E1399	JCM Folding Booster Base(Ht. Adj.)	\$2,786

Upper Extremity Support Tray (Angle Adjustable Hardware Included)

Tray Mount adjusts from 4-1/2"-7-1/2" with 1" tall cushion; Cut-out is 6"

☐	005A91	E0950	Light Grey	\$156
☐	005A92	E0950	Clear Poly	\$177

Shoe Holders & Foot Straps

☐	005A68	E0952, E0951	6" Long Shoe Holders - Pair	\$128
☐	005A69	E0952, E0951	8" Long Shoe Holders - Pair	\$128
☐	005A79	E0952	16" Foot Straps w/ Covers - Pair	\$38
☐	005A80	E0952	18" Foot Straps w/ Covers - Pair	\$38

Accessories

☐	005A74	K0105	IV Support Pole	\$156
☐	005A713	E2208	O2 Tank Holder	\$188
☐	005A714	K0105, E2208	IV/O2 Combination	\$263
☐	005LF5		Cup Holder/ Bottle/Feeding	\$48
☐	005LF7		Utility Hook (i.e. suction canister)	\$27
☐	005A462		Rain Cape	\$188



Front Mount Utility Platform



Cup Holder



Utility Hook



IV Support Pole



Upper Extremity Support Tray



Folding Booster Base



JCM Folding Booster Base



O2 Tank Holder




Battery and Vent/Utility Underframe Tray



Dovetail Vertical Vent Mount (LTV/Trilogy)



Rear Add On Accessory Platform



Child's Name _____

Age: _____

Diagnosis: _____

Height: _____

Weight: _____

1. Top of Head to Seat: _____

2. Top of Shoulder to Seat*: _____

3. Under Arm to Seat: _____

4. Trunk Depth: _____

5. Actual Sitting Depth*: _____

6. Foot Length: _____

7. Width of Head: _____

8. Shoulder Width: _____

9. Chest Width: _____

10. Hip Width: _____

11. Left Lower Leg Length*: _____

12. Right Lower Leg Length*: _____

* Critical Dimensions



Sunrise Medical (US) LLC
2842 Business Park Ave. · Fresno, CA 93727 · USA