

Delta® Pro Homecare Bed System

HCPCS: 15030N, 15004PN: E0295 | 15033NC, 15235N, 15005LPN:
E0297 | 15003PN: E0293,E0297, E0265, E0266



- Patent pending Universal Frame can configure with various component to create up to SIX unique bed configurations offering more homecare bed options without the need to store more beds!
- Assists with inventory constraints by offering 1 homecare bed system that can transform to suit multiple needs using color coded component options
- Product Safety & Performance Standards: Meets IEC 60601-2-52; 60601-1; 60601-1-2; 60601-1-11
- Safe working load capacity of 450 lbs
- Choose between semi-electric and full electric motor options to provide head, foot, and height adjustability of bed frame as needed (excludes 15003PN, manual crank homecare bed system)
- Ultra lightweight head and foot frames easy to set up by 1 individual with no links to connect and no tools required
- Square tube steel frame slat deck provides durable strength, no-sag mattress support, and is easy to clean
- HD Polyethylene molded bed ends are attractive, labeled for proper assembly, and easy to clean and maintain
- Bed End safety lock slides over hanger pins to prevent disengaging from bed ends
- UL approved motor is self contained for reduced weight and noise and can be easily installed and removed with patient in bed (excludes 15003PN)
- Reduce service calls with large glow in the dark hand pendant that is easy-to-use for adjusting head and foot deck sections; built in battery backup and LED indicators for power and lockout included (excludes 15003PN)
- Power failure protection available on dual style motors via a 9-volt battery (not included) located in the hand pendant which can return the frame to a horizontal position. Manual cranks included for all styles (excludes 15003PN, no electrical components means independent of electrical power availability)
- Frame includes mounted retractable full rail receiver brackets to eliminate cross braces ensuring proper placement and reducing equipment loss
- Half (15561HRAIL) and full (15561FRAIL) rails are available, sold seperately, and meet FDA entrapment guidelines as well as ASTM standards
-
-

| PRODUCTS | Description |
|----------|------------------------------------------|
| 15003PN | Delta Pro Bed - Manual |
| 15004PN | Delta Pro Bed-Standard Linear Semi Elect |
| 15005LPN | Delta Pro Low Linear Full Electric |
| 15030N | Delta Pro Bed -Standard Semi Electric |
| 15033NC | Delta Pro Bed - Standard Full Electric |
| 15235N | Delta Pro Bed -Low Full Electric |

| SPECIFICATIONS | Description |
|------------------------|------------------------------------------------------------------------------------------------------------------------------------|
| Battery | 15030N, 15033NC, 15235N: Backup available but not included |
| Bed Height | 15030N: 13.5" - 23" 15033NC: 13.5" - 23" 15235N: 10" - 23" 15003PN: 13.5" - 23" 15004PN: 13.5" - 23" 15005LPN: 10" - 23" |
| Bed Height w/3 Casters | 15030N, 15033NC, 15003PN and 15004PN: 28" - 38" 15235N and 15005LPN: 27" - 40" |
| Carton Shipping Weight | 15030N: 160 lbs 15033NC: 168 lbs 15235N: 171 lbs 15003PN: 166 lbs 15004PN: 173 lbs 15005LPN: 181 lbs |

Customer Service | 877.224.0946

Compressor Free Air Flow

Brown

4 Casters (2 locking and 2 unlocking)

ACCESSORIES**Description**

| | |
|------------|------------------------------------|
| 15550MOHC | Delta Pro Dual Motor w/pendant, EA |
| 15561FRAIL | Compet II Full Rail, 1/PR, PR |
| 15561HRAIL | Compet II Half Rail, 1/PR, PR |