

MEDMIZER®



## FLEXTILT CHAIR



For superior resident comfort and care.

# TODAY, WE ARE PROUD TO INTRODUCE THE *FLEXTILT CHAIR.*

At Med-Mizer, we have decades of experience and knowledge of the post acute care market. We also have a deep understanding of the challenges of being bedbound for the resident and caregiver.

We have built upon these past learnings and are excited to provide customers with an affordable solution, one that is Made in the USA, and from a name they have known to trust – Med-Mizer. Today, we are proud to introduce the FlexTilt chair.

## THE FLEXTILT WILL HELP YOU:

- Provide independence and mobility
- Create comfort with unlimited positioning options
- Improve skin/wound care with pressure relief
- Reduce caregiver injuries and efficiency (reduction in transfers, reduced time in repositioning, etc)
- Reduce cost with an affordable chair solution



NPUAP RECOMMENDS SHIFTING WEIGHT IN A CHAIR BOUND RESIDENT EVERY 15 MINUTES.<sup>1</sup>

## FLEXTILT CHAIR SPECIFICATIONS

<b>Frame</b> .....	Gauge steel 11ga - 16ga
<b>Degree of Tilt</b> .....	25°
<b>Degree of Recline</b> .....	30°
<b>Seat Width</b> .....	20", 22"*, 24"* *Longer lead times
<b>Seat Depth</b> .....	16.5"-18"
<b>Back Seat Height</b> .....	32"
<b>Seat Height</b> .....	19.5"-21.5"
<b>Weight Capacity</b> .....	350 LBS.
<b>Warranty</b> .....	3 year mechanical frame 1 year on components



### REFERENCES

1 | <http://www.npuap.org/wp-content/uploads/2014/03/NPUAP-F-tag-final-March-2014.pdf>

PROUDLY MADE IN THE U.S.A.

(877) 867-7365

**MED-MIZER.COM**



is at the forefront of innovating medical technology for post acute care. We are focused on partnering with you to help address healthcare challenges in the following key areas—




**FALLS  
PREVENTION**



**BARIATRIC  
CARE**



**MOBILITY**



**REHABILITATION**



**WOUND CARE**