


P2002 PRIME LONG TERM CARE BED



1. Head Deck- SP02-P2002-10003E

Head Deck Width Ext. Left- SP02-P2002-08E

Head Deck Width Ext. Right- SP02-P2002-07E
2. Foot Deck- SP02-P2002-10001E

Foot Deck Width Ext. Left- SP02-P2002-12E

Foot Deck Width Ext. Right- SP02-P2002-11E
3. Thigh Deck- Part no longer available.

Thigh Deck Width Ext.- SP02-P2002-09E
4. Hip Deck Width Extension - SP02-P2002-10E
5. Head/Foot/Thigh Deck Hdwe.- SP02-PPB2000-DB
6. Width Extension Guides- SP02-PPB2000-SMB
7. Width Extension Release Assy.- SP02-PPB2000-35E
8. Foot Deck Pivot Arm (pr)- SP02-PPB2000-27DB

Foot Deck Pivot Bolt (ea)- SP02-PPB2000-23DB
9. Wall Guard- SP02-PPB2000-32
10. Head/Foot Mattress Retainer- SP02-ML1002-10005E
11. Side Mattress Retainer- SP02-ML1002-10004E
12. Head Board Mtg. Bracket- SP02-PPB2000-1612E
13. Foot Board Mtg. Bracket- SP02-PPB1750-1712E
14. End Board Adj Mtg Brkt - SP02-PPB1750-31E

Hand Wheel Knob- SP02-PPB1750-43E
15. End Board Bolt Male (20)- SP02-UF16108051E

End Board Bolt Fem. (20)- SP02-UN90108091E
16. Total Locking Caster- SP02-UA18701500E
17. Direct Locking Caster- SP02-UA18701501E
18. Hand Pendant (8 button)- SP01-72298

Hand Pendant (10 button)- LA-IPROXX2TREND
19. Motor Control Box- SP01-69159
20. Motor Control Box Bypass Plug- SP01-PCB900-47E
21. Motor Control Box Plug Cover- SP01-PCB601-01-53E
22. Power Cord- SP01-CA65738
23. "Y" Cable - SP01-CA51336
24. "Y" Cable Mounting Bracket- SP02-PPB2000-36E
25. Scale (PBAWS & PASCB)- SP02-PPBASC [1]
26. Hi/Lo Motor (2 used)- SP01-79323

Hi/Lo Pivot Bolt & Hdw. (2") (ea)- SP02-ML10021000701
27. Foot Motor- SP01-81069
28. Head Motor- SP01-81077
29. Side/Assist Rails- SP02-PPBASR [1]

Side/Assist Rail Pad- SP02-PPBASRP
30. Side/Assist Rail Mounting Bracket- SP02-PPBRSRB

Side/Assist Rail Rel. Handle- SP02-PPB1750ARRH

Side/Assist Rail Liner- SP02-1750-13-03E
31. Head/Foot Board Plug (pr)- SP02-PPB2000-17-03
32. Frame Tube End Insert- SP02-PPB2000-42E
33. Lt Side Head/Rt Side Brake Pedal- SP02-PPB2000-10019

Rt Side Head/Lt Side Brake Pedal- SP02-PPB2000-10018

Plastic Pedal Washer- SP02-PF43131900E

Brake Pedal "E" Clip- SP02-UE31091001E

Brake Pedal Sm. Flange Washer- SP02-UB10122501E

Brake Pedal Lg. Flange Washer- SP02-UB10122500E

Brake Pedal Ext. Bar- SP02-ML1002-10012E
34. Cross Bar Mounting Plate- SP02-PPB2000-21B

Cross Bar Mtg. Bushing- SP02-PPB2000-22B
35. Cross Brake Locking Assy.- SP02-PPB2000-24B
36. Length Extension Bolt- SP01-KBCZ31C50
- n/s Hardware Pack- SP02-ML1000-10030E
- n/s Embedded Staff Control- [2]

[1] Accessory, not a part.
[2] No longer available