

STAR 6

4256

The NOVA STAR shines with more features, options, colors and style. Its *Quick-Fit* adjustable design gives the user a perfect fit in just seconds; also making it so easy to pack up for travel, adventure and those odd spaces.



18.5", 19.5",
20.5" & 21.5"
Adjustable
Seat Height

COMES STANDARD

- Locking Hand Brakes
- Quick-Fit* Push-button Adjustable
- Padded Seat w/ Under-Seat Pouch
- Built-in Foot Pedals
- Oval-shaped Tubing
- Extended Folding Handle
- 6" Wheels

OPTIONS

- Cup Holders (CH-1100 & CH-3100R)
- Cane Holder (CH-4100R)
- Designer Mobility Bags
- Back and Seat Covers
- Travel Bag

Colors: Black, Blue, Diamond Sky Blue, Purple, Red

WEIGHT CAPACITY: 300 lb

SEAT HEIGHT: 18.5", 19.5", 20.5" & 21.5"

SEAT DIMENSIONS: 14.25"w x 14.25"d

HANDLE HEIGHT ADJUSTMENT: 30.75" - 38.25"

WIDTH BETWEEN HANDLES: 18.5"

OVERALL DIMENSIONS: 24.5"w x 28.75"d

APPROX. USER HEIGHT: 5'1" - 6'2"

POUCH DIMENSIONS: 14"w x 9"d x 6.25"h

DIMENSIONS WHEN FOLDED: 24.5"w x 45"d x 12"h

WEIGHT: 14.75 lb



4258

STAR 8

The NOVA STAR shines with more features, options, colors and style. Its *Quick-Fit* adjustable design gives the user a perfect fit in just seconds; also making it so easy to pack up for travel, adventure and those odd spaces.



19.5", 20.5"
& 21.5"
Adjustable
Seat Height

COMES STANDARD

- Locking Hand Brakes
- Quick-Fit* Push-button Adjustable
- Padded Seat w/ Under-Seat Pouch
- Built-in Foot Pedals
- Oval-shaped Tubing
- Extended Folding Handle
- 8" Wheels

OPTIONS

- Cup Holders (CH-1100 & CH-3100R)
- Cane Holder (CH-4100R)
- Designer Mobility Bags
- Back and Seat Covers
- Travel Bag

Colors: Black, Blue, Diamond Sky Blue, Purple, Red

WEIGHT CAPACITY: 300 lb

SEAT HEIGHT: 19.5", 20.5" & 21.5"

SEAT DIMENSIONS: 14.25"w x 14.25"d

HANDLE HEIGHT ADJUSTMENT: 32" - 38.5"

WIDTH BETWEEN HANDLES: 18.5"

OVERALL DIMENSIONS: 25"w x 27.25"d

APPROX. USER HEIGHT: 5'3" - 6'2"

POUCH DIMENSIONS: 14.5"w x 8.75"d x 6.25"h

DIMENSIONS WHEN FOLDED: 25.5"w x 40.5"d x 15.5"h

WEIGHT: 16.5 lb

