WHEEL-A-BOUT WHEEL-A-BOUT





- Closure strap locks ramp panels together
- Easy-to-handle & set up
- Durable welded construction
- Full platform provides excellent stability
- Anti-slip, high traction surface
- Accommodates wheelchairs and scooters with various wheel configurations
- Stores upright in garage
- Safety DVD and steel security pins included

Effortless Wheel Design



Quick Release Pins used to secure Latch Straps

WAB1030 & WAB1230 Feature 2 3/4" Long Top Hook

Solid / Stable Surface



| Item Number | Top Hook Length | Length | Width | Folded Size | Weight | Capacity* |
|--|-----------------|--------|-------|---|---------|-----------|
| WAB1030 | 2 3/4" | 10' | 30" | 66 ¹ / ₄ " × 36 ³ / ₄ " × 16" | 84 lbs. | 600 lbs. |
| WAB1230 | 2 3/4" | 12' | 30" | 78½" x 36¾" x 16" | 98 lbs. | 600 lbs. |
| WABU1030 | 9 3/4" | 10' | 30" | 63½" x 36¾" x 16" | 83 lbs. | 600 lbs. |
| WABU1230 | 9 3/4" | 12' | 30" | 75¼" x 36¾" x 16" | 97 lbs. | 600 lbs. |
| Not sure which ramp length you need? Refer to pages 2 - 3 to help you. | | | | | | |

^{*} Maximum ramp capacity: 300 lbs. (one axle) – 600 lbs. (two axles).



