

VivaLIFT![®]
POWER RECLINERS
By Pride[®] Mobility Products Corp.

THE MOST
COMFORTABLE
RECLINER
IN THE WORLD!

THE 4-MOTOR RECLINER

WITH
COUNTLESS
POSITION
COMBINATIONS



Shown in Badlands
Walnut Fabric

ELEGANCE

A truly elegant, sophisticated recliner is more than just physically alluring: It's the precise combination of function, style and comfort. The Elegance Collection from **VivaLift!**[®] Power Recliners is all of these and more, providing a cozy place to relax with your family. The soft fabric and contrast stitching makes this recliner a perfect addition to your home.



Pride Mobility[®] (US) 800-800-8586 | (CAN) 888-570-1113
pridemobility.com | prideliftchairs.com |    

ELEGANCE

COLLECTION	ELEGANCE	ELEGANCE
MODEL	PLR-975M	PLR-975L
POSITION	Infinite Lay Flat Relaxed Reclined Upward Lift	Infinite Lay Flat Relaxed Reclined Upward Lift
SIZE	Medium	Large
WEIGHT CAPACITY	375 lbs.	375 lbs.
SEAT-TO-FLOOR	20"	20.5"
SEAT DEPTH	21"	22"
SEAT WIDTH	21"	22"
TOP OF BACK TO SEAT	27"	28"
FULLY PADDED CHAISE	Yes	Yes
FOOTREST EXTENSION	Standard	Standard
SEATING TYPE	Total Comfort	Total Comfort
CHAIR HEIGHT	42"	43"
DISTANCE FROM WALL	14"	16"
OVERALL WIDTH	35.5"	37.5"
CHAIR WEIGHT	152 lbs.	157 lbs.
WARRANTY	7-Yr. Pro-rated	7-Yr. Pro-rated

FUNCTIONS AND FEATURES:

-  POWER HEADREST
-  POWER LUMBAR
-  USB REMOTE
-  LITHIUM BATTERY BACKUP
-  FOOTREST EXTENSION
-  POWER BACKREST

POSITION GUIDE:

-  RELAXED
-  RECLINED
-  UPWARD LIFT
-  INFINITE LAY FLAT

STANDARD FABRIC: BADLANDS



MUSHROOM

WALNUT

STEEL

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

 = Indicates shortest available lead time

Recommended HCPCS: E0627-Seat lift mechanism, electric, any type DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. Payer coding, coverage and bundling guidelines may apply.



Pride Mobility® (US) 800-800-8586 | (CAN) 888-570-1113

pridemobility.com | prideliftchairs.com |    