PET CARRIER & CUSHION BLACK

DON'T THINK TWICE,
JUST BRING IT ALONG!



movinglife[®]

ATTO Accessories are created specifically to enhance, support and work with the ATTO mobility scootesr in all user modes; while driving, while wheeling along as trolley and when folded for storage or transportation.



SEATING & SUPPORT



ATTO /ATTO SPORT Armrests



Seat cushion



Seat cushion leather



Luxury Swivel Seat Leather



Belt & Support Strap



Luxury Orthopedic Backrest

STORAGE



Essentials Pouch



Front & rear basket



Carryall & Cushion



Carryall & Cushion leather



Pet Carrier & Cushion Black



Flight Kit



Oxygen Tank Holder and Case

TECHNICAL



Cane / Crutches Holders



Mobile Phone Holder



Rear Light Lifting Handle



ATTO Left Hand Throttle Kit



Standard Flight Battery




Launch – H2 2023 Long distance battery – not flight

approved



Cable kit

How to Use PET CARRIER & CUSHION BLACK



Slip the cushion over the seat (use handles to pull if needed)




02

Simply pull hook on elastic strap to fold under and connect to ATTO chassis

NOTE!

Take off before folding ATTO



03

One-hand "pull to open" secure funnel access.

Quick "no zip" flap, to ensure your pet stays inside when closing.



04

Attach to ATTO Seat Cushion dedicated upper zip