

ziggo Pro

PRODUCT FEATURES



SPECIFICATIONS	ZREC1200	ZREC1400	ZREC1600
	VE CARRIE SPECIAL SA	The Control of the Co	
FRAME COLORS	ORANGE	RED	BLUE
SEAT WIDTH	12"	14"	16"
SEAT DEPTH	12"	14"	16"
OVERALL OPEN WIDTH	20"	22"	24"
OVERALL FOLDED WIDTH	12"	12"	12"
BACK HEIGHT INCL HEADREST	25"	25"	31"
BACK HEIGHT W.O HEADREST	17.5"	17.5"	19.5"
OVERALL HEIGHT	40.75"	40.75"	49.5"
SEAT TO FLOOR HEIGHT	15"/16"	15"/ 16"	19.5"/ 20.5"
OVERALL LENGTH W/WRIGGINGS	38.5"	38.5"	52.75"
BACKREST RECLINES	90°- 180°	90° - 180°	90°- 180°
NET WEIGHT	44 LBS	44 LBS	56.75 LBS
CARTON SHIPPING WEIGHT	50 LBS	50 LBS	62 LBS
WEIGHT CAPACITY	150 LBS	200 LBS	300 LBS







Colored Steel Frame, Lightweight & Foldable Design



Powered by Circle Specialty's 25+ years of experience in revolutionary mobility and rehab equipment tailored exclusively to the needs of kids and teens, you have every good reason to make Ziggo Pro a part of your family.