

ziggo



LET THE ziggo TAKE YOU PLACES

INDOORS AND OUTDOORS



THE LIGHTWEIGHT, FOLDING
MANUAL WHEELCHAIR
TO TAKE YOU PLACES



AVAILABLE IN 4 SIZES:

12", 14", 16" & 18"



DESIGNED FOR EASY
TRANSPORT & STORAGE



THE ziggo IS A LIGHTWEIGHT MANUAL WHEELCHAIR

designed for kids, teens and adults **on-the-go**. Loved by users and caregivers alike, the ziggo is practical, durable, adjustable and cost effective, allowing for maximum maneuverability and easy transport, be it indoors and outdoors. Meeting the standards of therapists and clinicians, the ziggo provides comfortable support, functionality and positioning options.

Experience the little things in life with great ease.

FEATURES

Colored steel frame

Folds easily for compact storage and transport

Comfortable nylon upholstery with back pocket

Flip-back padded desk length arms

Fold-down backrest on 16" and 18" Ziggo model.

Pelvic belt

Quick-release treaded PU rear wheels

Solid PVC front wheels

Push-to-lock brakes

Height adjustable swing-away footrests with pastic footplates

Calf strap for extra support and comfort

Variety of optional accessories to meet your individual needs

3-year limited warranty



SPECIFICATIONS

	ZG1200	ZG1400	ZG1600	ZG1800
Seat Width	12"	14"	16"	18"
Seat Depth	12"	16"	16"	16"
Back Height	14"	16"	16"	19"
Seat Height Front	18"	18.5"	18.5"	19"
Seat Height Back	17"	17.5"	17.5"	18"
Overall Open Width	18.5"	21"	23"	25"
Overall Length w/Riggings	36"	41"	42"	42"
Weight w/o Riggings	28 lbs	36.5 lbs	37 lbs	39 lbs
Weight Capacity	150 lbs	150 lbs	250 lbs	250 lbs

Wheelchair accessories play a vital role in positioning and comfort and enhance the functionality, support and safety of your wheelchair. They can be configured based on your needs and preferences.

Installation is simple, giving you the ease to add or remove as required.



Height, length and angle adjustable swing-away elevating legrests with padded calf pads.

ZGELR01 for use with ZG1200 & ZG1400

ZGELR02 for use with ZG1600 & ZG1800

ELEVATING LEGRESTS

Padded trunk harness offers anterior torso support and improved posture alignment. Easily adjustable for growth and change of positioning requirements.

ZGTH01 (use with all Ziggo products)



TRUNK HARNESS



Increases wheelchair stability and prevents wheelchair from tipping over backwards.

ZGAT01 (use with all Ziggo products)

ANTI-TIPPERS

Height, depth and angle adjustable headrest provides posterior head and neck support.

ZG8100 (use with all Ziggo products)



HEADREST



Lightweight foam cushion improves comfort and sitting tolerance.

ZGCS12 for use with ZG1200

ZGCS14 for use with ZG1400

ZGCS16 for use with ZG1600

ZGCS18 for use with ZG1800

SEAT CUSHION