

# Klip



**CHOICES  
PROVIDE YOU  
WITH THE  
BEST OPTIONS**



The **Klip** is a lightweight posterior gait trainer for kids, teens and young adults with special needs who require support and assistance with their day-to-day walking.

Placed behind the user, the Klip is open in the walking position and allows freedom of movement and promotes walking in an upright position. It is recommended for those with Cerebral Palsy and other neurological conditions.

The **Klip** can also be used in the anterior position for all sizes other than the XS size).

It is simple and safe to use, easy to fold and can be used as a basic posterior gait trainer or as a fully supported gait trainer.

**SPECIFICATIONS**

	Width inside Handles	Handle Height	Seat to Floor Height	Outside Dimension (Open at Base)	Product Weight	Weight Capacity
<b>2 wheeled</b>						
KP210	13.5"	16"-195"		23" x 23"	7 lbs.	75 lbs.
KP220R	14"	21"-27"		24" x 31"	9 lbs.	85 lbs.
KP230	16"	24"-32"		23" x 35"	11 lbs.	150 lbs.
<b>4 wheeled</b>						
KP410	13.5"	16"-195"		23" x 23"	9 lbs.	75 lbs.
KP420R	14"	21"-27"		24" x 31"	11 lbs.	85 lbs.
KP430	16"	24"-32"		23" x 35"	13 lbs.	150 lbs.
KP440B	17"	30"-37"		26" x 37"	13 lbs.	190 lbs.
<b>4 wheeled with Seat</b>						
KP510	13.5:	16"-195"	9.50" - 13.50"	23" x 23"	11 lbs.	75 lbs.
KP520R	14"	21"-27"	12" -18"	24" x 31"	13 lbs.	85 lbs.
KP530	16"	24"-32"	14" -21"	23" x 35"	17 lbs.	150 lbs.
KP540B	17"	30"-37"	17" -24"	26" x 37"	18 lbs.	190 lbs.
<b>4 wheeled with 8" wheels</b>						
KP418						
KP428R	14"	23"-30"		24" x 33"	15 lbs.	85 lbs.
KP438	16"	28"X35"		25" x 37"	17 lbs.	150 lbs.
KP448B	17"	33"X41"		26" x 40"	18 lbs.	190 lbs.
<b>4 wheeled with seat and 8" wheels</b>						
KP518						
KP528R	14"	23"-30"	15" -21"	24" x 33"	17 lbs.	85 lbs.
KP538	16"	28"X35"	17" -24"	25" x 37"	21 lbs.	150 lbs.
KP548B	17"	33"X41"	20"X28"	26" x 40"	21 lbs.	190 lbs.



### 2 Wheeled

The 2-Wheeled Klip is an advantage for users with significant balance impairments since its two rear crutch tips provides extra resistance and is slower moving than a 4 Wheeled walker.

*It is available in three sizes for toddlers, kids and young adults.*



### 4 Wheeled

The 4-Wheeled Klip with front swivel wheels makes it easy to turn and provides maximum maneuverability.

*It is available in four sizes to accommodate toddlers, children, teens and adults.*



### Outdoors

The 4-Wheeled Klip with 8" swivel wheels is easy to navigate on all indoor and outdoor surfaces. It is excellent for an active lifestyle.

*It is available in four sizes to accommodate toddlers, children, teens and adults.*



### Flip-up seat

The 4-Wheeled Klip with flip-up seat is a great convenience for users who may lack the endurance to walk for long periods of time. The seat flips up for walking and standing and folds down for sitting and resting.

*It is available in four sizes to accommodate toddlers, children, teens and adults.*



### Outdoors/Seat

Excellent on many fronts

Simple to use and maneuver, 8" wheels for smoother rolling on both indoor and outdoor surfaces and a seat that folds down for sitting and resting.

*It is available in four sizes to accommodate toddlers, children, teens and adults.*

**Height adjustable to ensure proper fit and accommodate future growth**  
**Foot-Activated swivel wheel lock with an audible "click" for easier and greater control.\* On/Off one-directional rear wheels prevents walker from rolling backwards. Simple to disengage when walking skills improve\*\***  
**Folds easily for compact storage and transport**

*\*Not available on the 2 wheeled and outdoor models*

*\*\*Not available on the 2 wheeled models*

## ACCESSORIES

The **Klip** is a basic gait trainer with a choice of accessories designed to improve support, posture and mobility and offer individualized assistive support.

They can be easily added or removed as their walking skills and needs evolve.



**FOREARM PLATFORM WITH HANDGRIP**

Fully adjustable forearm platforms are height and depth adjustable and the handgrips extend forward or back, supinate, pronate, rotate and are wrist-angle adjustable.



Extends



Rotates



Wrist  
Angles



Wrist  
Rotation

Width and height adjustable with a support belt.



**PELVIC STABILIZER**



**SLING SEAT**

Adjustable cushioned seat provides partial weight bearing and properly positions user's pelvis.