MATTRESS REDISTRIBUTION with Memory Foam







REDISTRIBUTION FOAM MATTRESS - STAGE II MAXIMUM PREVENTION FOR PRESSURE ULCERS



HIGH DENSITY Durable and Flexible



ANTI SHEAR SURFACE





WEIGHT CAPACITY: 300 Lbs.



HEEL SLOPE for maximum wound prevention at the heel



Beveled Perimeter Edge allows for safe and easy transfer in and out of bed

HOSPITAL GRADE COVER

- 70 Denier Nylon with PVC backing
- Detachable cover with a protected zipper for sanitary benefits and to prevent infections
- Easy to maintain
- Waterproof
- Anti Microbial
- Anti-Flammable Available with a Fire Sock
- Flame Retardant CFR1632 1633
- Conforms to CAL TB 117



MADE FROM QUALITY EUROPEAN FOAM

MODEL # PMX8036 80" x 36" x 6"

www.medacureinc.com