4-Wheel Full Suspension Electric Scooter

- > Brand new sporty streamline design
- > Full suspension system with shock absorbers delivers a smooth and comfortable ride
- Modern LED lighting system delivers a brilliant light with less power consumption
- Digital LCD dashboard display: time display, battery gauge, temperature, maintenance messages, speed display and odometer
- > Extra storage space beneath tiller
- Automatic safety turning half speed control by sensor





LED lighting system



LED lighting system



Digital LCD dashboard display




More space for storage

Technical Specs.

Weight Capacity	159kg/ 350lb					
Max. Speed	12kph/ 7.5mph; 15kph/ 9.5mph					
Range	50km/ 32mi.					
Turn Radius	167cm/ 65.7"					
Ground Clearance	9cm/ 3.5"					
Motor	DC 24V, 400W(Rating)/ 1750W(Max)					
Controller	120A					
Battery	12V/ 75AH x 2pcs					
Charger	5A off-board					
Gradient	up to 10°					
Front Wheel	12"					
Rear Wheel	12"					
Brake	intelligent, regenerative, electromagnetic brakes					

Dimensions







seat dimensions [cm/in]										
D	E		F	G	G		Н		I	
46/18	40,	/16	46/18	63.5/25		43-48/17-19		97/38		
overall dimensions [cm/in] weig						ht [kg/lb]				
А	В	С	base	battery/ea		seat	ttl w/batt.		capacity	
130/51	69/27	125/49	69/152	25/55	1	8/40	118/260		159/350	