

Prime Plus Bed Model P2002

HCPCS:




FIGURE A



FIGURE B



FIGURE C



FIGURE D

- Quiet, smooth operation/Easy-to-use Hand Pendant electrically controls Head/Foot/Contour positioning and overall height of the bed frame for maximum comfort and versatility for the patient and caregiver
- Clip-on pendant mounts on bed rail, within easy reach to operate all adjustments
- Variable Comfort Chair and Cardiac Chair (with Trendelenburg/Reverse Trendelenburg staff control option) adjustability helps provide added comfort and relaxation positioning for the resident and facilitates caregiver accessibility
- Cardiac Chair positioning is beneficial for residents with heart and respiratory problems that find the semi-upright position brings comfort and relief by elevating the torso from the waist up, which provides relief to the lungs, increases circulation
- Trendelenburg/Reverse Trendelenburg capable positioning with Staff Control Option is ideal for rehab and sub acute care, provides therapeutic benefits that can aid in circulation and ease discomfort
- Edema foot ratchet elevates the foot end of the bed
- 4 heavy duty Euro style roll at any height steer & locking casters provide easy maneuverability and secure lock positioning
- Adjustable wall saver minimizes distance between wall to bed frame to prevent damage with/without optional Trapeze
- Adjustable mattress retainer securely supports mattress position for 80", 84", 88" length
- Built in push button width expansion from 39" to 54" width
- All Steel frame and no sag slat deck construction provides superior strength, durability, mattress support and is easy to clean
- Optional oversized Bed Ends with bumper bar available in choice of 5 attractive wood grain laminate colors will compliment most any décor
- One button drop down padded Assist Rails are standard for patient comfort and safety
- Optional Embedded Staff Control at the foot end of the bed with Trendelenburg/Reverse Trendelenburg positioning and lock out feature provides easy access for position adjustments
- Optional extra padded assist bars, trapeze and integral weigh scale (P2002-WS)
- Optional Battery Backup for emergency use when loss of power occurs
- 1000 lbs weight capacity*
- 80", 84", 88" sleep surface
- 39" width w/integral expansion to 54"
- 19" low height 32" high working height
- Intertek ETL Classified
- Certified to CAN/CSA Std C22 2 No 60601-1 & 60601-2-52
- Conforms to ANSI/AAMI Std E60601-1, IEC Std 60601-1, IEC STD 60601-2-52
- International Protection (IP) 54
- CE Safety Certification
- Warranty: 15 years on frame lifetime on welds 3 years on electronics
- *Item #P2002-WS includes integral weigh scale option installed onto the bed frame

PRODUCTS

Description

P2002

PrimePlus Bariatric Bed/1000 lb. Wgt cap

SPECIFICATIONS	Description
Base Width	37"
Built-in Width Expansion	39", 42", 48", 54"
Sleep Surface	80" - 88"
Super Low Height	19"
Trendelburg/Reverse Trendelenburg Capable	With Staff Control Option
Warranty	15 Year on Frame, 3 Year on Electric and Lifetime on All Welds
Weight Capacity	1000 lbs *Incl. resident, mattress, bed ends, assist bars, all access, personal items
Working Height	32"
ACCESSORIES	Description
PBABE-FM	BARI BE,FMAHOG,F/1750/52,2000/02,EA, EA
PBABE-B	BARI BE,OAK,F/1750/52,2000/02,EA, EA
PBABE-A	BARI BE,CHERRY,F/1750/52,2000/02,EA, EA