



Expandable Bariatric Bed




The Infinity bed is a fully electric bed designed with safety in mind for the diverse needs of the residents. With integrated width and length expandability, 600 lb capacity, low position, and many more standard features, the Infinity bed will become every facility's preferred bed by residents and staff alike.



**INTEGRATED
EXPANSION**



**COMFORT CHAIR
POSITION**



**MOBILITY AT
ANY HEIGHT**

SPECS

Model	54000
Bed Width Range	Integrated 36" – 42" – 48"
Bed Length Range	Integrated 80" – 84"
Bed Height Range	7" – 30"
Weight Capacity	600 lbs.
Warranty	15 yrs. Welds/ 3 yrs. Electric/ 1 yr Components

FEATURES

- Integrated width/length expansion
- Comfort chair position
- Auto contour
- Battery backup
- Mobility at any height
- Central locking system
- Quick release head & footboards
- Underbed lighting
- Bumper guard
- Safety stop
- Heel lift extension




Expandable Bariatric Bed

✉ info@emeraldsupply.net | 🌐 emeraldsupply.net



Bumper guard




Integrated length expansion



One touch central lock




Integrated width expansion



Auto Contour



Comfort chair position



RAILS



Q Assist Bar



Pivot Assist Bar



Soft Style Side Rails

BARS & RAILS	ITEM #
Soft Style Side Rails	52210
Q Assist Bar	52211
Pivot Assist Bar	52212

HEAD & FOOTBOARDS



Solar Oak



Mahogany



Amber Cherry

BOARD COLORS	ITEM #
Solar Oak	54020
Mahogany	54021
Amber Cherry	54022