

Be Independent.™



NEW

EZ STAND-N-GO™

BECAUSE STANDING FROM YOUR FAVORITE COUCH OR CHAIR CAN BE EASY



ERGONOMIC HANDLES
FOR EXTRA COMFORT



WATCH A DEMO



PRODUCT #2200





Universal Fit - Even When Reclining!



Fits Under Cushion of Couch or Chair



Fits on Top of Couch, Chair or Recliner Cushion

FEATURES & BENEFITS

- Stand safely from any couch, chair, or recliner
- Dual support cushion handles for easy standing
- Universal fit: on top of, or between chair cushions
- Fully adjustable in height and depth

SPECIFICATIONS

Chair Type: Fully adjustable to fit most couches, chairs, and recliners

Handle Height From Floor: Adjustable from 22" - 32"

Width Between Handles: Adjustable from 19" - 26"

Base Dimensions: Adjustable Depth 23" - 26" and Adjustable Width 19" - 26"

Tube Diameter: 1"

Product Dimensions: 23" - 26" x 19" - 26" x 22" - 32" 8.5 lbs

Package Dimensions: 21" x 12.5" x 5.5" 8.5 lbs

Case of 5 Dimensions: 25.5" x 21.25" x 12.5" 42.5 lbs

Made Of: Powder Coated Steel

Weight Capacity: 300 lbs