



Invacare® Top End® Schulte™ 7000 Series Basketball Wheelchair

The **Invacare® Top End® Schulte™ 7000 Series Basketball Wheelchair** is faster, lighter and stronger.* In 2011, it hit the market with new material, enhanced headtube design and axle securing system, reducing overall weight by two pounds and improving performance. Each Schulte 7000 Series Basketball Wheelchair is custom built to each individual athlete's specifications and is overseen by Paul Schulte who is the design engineer for Invacare Top End and three-time Paralympian as a member of the USA National Men's Wheelchair Basketball Team.


*As compared to previous Invacare® Top End® model wheelchairs.

Features

Custom-built fixed frame in 7005 aluminum
Adjustable height fixed swivel anti-tip 59mm fifth wheel
Fixed back with fixed or adjustable height
Adjustable tension seat and back upholstery
Fixed camber: 15°, 18° or 20°
High performance Sun® wheels: 24", 25", 26" or 700C
Four bend offensive wing

High pressure clincher tires
Precision fork caster system
72mm casters on front, 59mm casters on rear
Powder coat paint
Rear frame protector
250 lb. weight capacity with optional heavy duty upgrade for 350 lb. capacity

Model #: PS7



1. Double swivel anti-tip/ fifth wheel option
2. Precision fork caster system
3. Optional welded fixed tube sideguard
4. Adjustable center-of-gravity option
5. Adjustable height fixed swivel anti-tip
6. 7005 Aluminum heat treating process



Invacare Corporation
www.invacare.com

USA

One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada

570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

Top End

www.topendwheelchair.com
Sports and Recreation
Products
4501 63rd Circle North
Pinellas Park, Florida
33781-5914
(800) 532-8677
Fax: (727) 522-1007

Options

Adjustable Center-of-Gravity
Spinergy® SLX and XSLX wheels
Sun® Equalizer™ wheels
12 Tab pushrims
Tapered seat width
Welded seat rigidizer bars with solid seat
Double swivel anti-tip/fifth wheel
Click straps

Recessed platform footrest
Adjustable angle footrest
Clamp-on fixed aluminum side guards
Fixed side guards with welded tube
Fixed aluminum or carbon fiber side guards
Front frame protectors
Alignment gauge

©2018
Invacare Corporation.
All rights reserved.
Trademarks are identified
by the symbols ™ and ®.
All trademarks are owned
by or licensed to Invacare
Corporation unless
otherwise noted. Spinergy
is a registered trademark
of Spinergy Inc. Sun
and Equalizer are
registered trademarks
of Sun Components.
Specifications are subject
to change without
notification.
Form No. 18-179
180737

TOP END®