

Be Independent.™



# SECURITY POLE™

BECAUSE YOUR SECURITY POLE CAN GO ANYWHERE IN THE HOUSE



PRODUCT #1150



WATCH A DEMO





## FEATURES & BENEFITS

- **Comfortable Cushion Grip** – Easy-to-sanitize closed cell foam provides non-slip surface.
- **Adjustable** – Fits ceiling heights 7–10 ft.
- **Portable** – Main pole easily separates into two 53" pieces making it easy to transport.
- **Adaptable** – No need for wall mounts; can be installed in virtually any area of a room exactly where it is needed.

## SPECIFICATIONS

*Materials:* Rust resistant zinc plated steel with durable powder coat finish

*Fits Ceiling Heights:* Range of 7-10 ft.

*Weight Capacity:* 300 lbs

*Weight of Product:* 20 lbs

*Package Dimensions:* 54.00" L x 12.50" W x 3.25" H; 22 lbs

*Case of 2 Dimensions:* 54.00" L x 12.50" W x 6.5" H; 45 lbs

## COMES IN TWO COLORS

■ BLACK #1150-B

□ WHITE #1150-W

## EASY ASSEMBLY (No Additional Tools Required)



Attach Top and Bottom Pieces Using Included Pins



Slide the Two 53" Pieces Together



Use Included Wrench to Tension Mount Between Floor and Ceiling