



Wheelchair

The iCruise standard wheelchair is designed to meet a variety of patient requirements. It encompasses an array of features, including extra padding, comfortable vinyl upholstery, and elevating legrests.



FEATURES

- Slide tube, powder-coated silver vein steel frame
- For the wheelchairs with widths of 16"-20", dual axle provides easy transition of seat height to hemi-level
- Easy-to-clean, comfortable vinyl upholstery resists mildew
- Padded, upholstered removable desk armrests
- Includes elevating legrests with padded calf pads
- Elevating legrests are equipped with a pneumatic piston for infinite elevating positions
- Easy push-to-lock wheel locks
- Precision sealed wheel bearings in front and rear ensure long lasting performance and reliability
- Composite handrims with 24" solid urethane tires
- 8" x 1" solid PVC front casters



SPECS

Item #	Seat Width (in)	Seat Depth (in)	Back Height (in)	Floor to seat Height (in)	Rear Wheel	Overall Length (in)	Overall Width (in)	Overall Height (in)	Overall Width Folded (in)	Weight Capacity
11916	16	16	18	19.5	24"	45.5	24	36	12	300 lb
11918	18	16	18	19.5	24"	45.5	26	36	12	300 lb
11920	20	16	18	19.5	24"	45.5	28	36	12	300 lb
11922	22	18	18	21	24"	48	30	36	13	350 lb
11924	24	18	18	21	24"	48	32	36	13	450 lb